

PORTABLE XRD ANALYZER - X-RAY DIFFRACTOMETER - SHINE RESEARCH EDITION



The Portable XRD Analyzer – X-ray Diffractometer – SHINE Research Edition is integrating XRD, XRF, and computer software technologies. It enables the analysis of samples through X-ray diffraction, providing insights into material composition and internal structure, including crystal phase identification and trace elements.

With simple sample preparation and high precision measurements, it offers a comprehensive understanding of crystal integrity, making it a primary method for studying phase and crystal structure in materials.

Utilizing transmission geometry technology, SHINE directs X-ray beams through the sample location. A micro-focus X-ray tube emits beams through a collimator, irradiating powder particles in the sample pool. A unique vibration system ensures increased exposure of crystal powder particles to the X-ray beam, adhering to Bragg's law. The fixed sample size in the pool maintains resolution stability, crucial for low-density materials like drug samples, where higher resolution is achieved compared to instruments utilizing reflective geometry.

APPLICATIONS

The Portable XRD Analyzer – X-ray Diffractometer – SHINE Research Edition supports bio-pharmaceutical research, environmental monitoring, geology and prospecting, metallic compound chemistry, catalyst development, and teaching laboratories.

Portable XRD Analyzer - Bio-pharmaceutical

Bio-pharmaceutical

Portable XRD Analyzer - Environmental Monitoring

Environmental Monitoring

Portable XRD Analyzer - Geology and Prospecting

Geology and Prospecting

Portable XRD Analyzer - Research and Education

Research and Education

Portable XRD Analyzer - Metallic Compound

Metallic Compound

Portable XRD Analyzer - Chemistry and Catalysts

Chemistry and Catalysts

FEATURES

Portable

The Portable XRD Analyzer – X-ray Diffractometer – SHINE Research Edition features an all-in-one waterproof and dustproof design, eliminating the need for mechanical moving parts. It is lightweight, compact, and easy to carry, facilitating seamless scientific research in both laboratory and field settings.

Integration

XRD is seamlessly integrated with XRF technology, allowing for simultaneous collection of X-ray photon data from both methods during each test. This comprehensive approach provides detailed information on material composition, phase, and structure, thereby improving the accuracy of test results.

Easy-Operation

The Portable XRD Analyzer – X-ray Diffractometer – SHINE Research Edition is user-friendly, requiring no calibration and offering automatic detection. It provides real-time display of full spectrum diffraction, automatic phase analysis, and generates quantitative results with just one key operation.

Sample preparation is simple

For a single analysis, only about 20 mg of sample is required to achieve high-quality test results. Sample preparation, including crushing, filtering, and sampling, can be completed in just 3 minutes without the need for molding, pressing, or scraping for uniformity.

Data transmission

The Portable XRD Analyzer – X-ray Diffractometer – SHINE Research Edition achieves high-speed connection with a laptop computer via USB, Bluetooth, and WIFI, enabling real-time control and phase analysis of SHINE.

Safety

The Portable XRD Analyzer – X-ray Diffractometer – SHINE Research Edition undergoes professional multiple protective radiation treatments to ensure there is no radiation leakage in any direction during measurement.

THEORY & METHOD

X-ray diffraction (XRD) identifies the crystalline structure of a material, which reveals what phases or compounds it contains. Because every crystalline substance has a unique atomic arrangement, XRD acts like a fingerprint reader for solids.

Bragg's Law and Diffraction


The analyzer directs a beam of X-rays at the sample. The atoms in a crystal are arranged in regular planes, and these planes reflect the X-rays. When the path difference between rays reflected from adjacent planes equals a whole number of wavelengths, the rays reinforce and produce a strong diffracted beam. This condition is described by Bragg's law, where n times the wavelength equals $2d$ times the sine of θ , and d is the spacing between crystal planes.

From Diffraction Pattern to Identity

The instrument measures the angles and intensities of the diffracted beams to build a diffraction pattern. Each crystalline phase produces a characteristic set of peaks, so comparing the pattern against reference databases identifies the phases present and can reveal crystallite size and strain. This makes XRD a core tool for material identification and quality control.

TECHNICAL SPECIFICATIONS

XRD resolution	0.2°@2θ FWHM
Scope of XRD	5-55°2θ
The detector	2000 X 256 pixel, 2d, level 3 Pelter refrigeration CCD
X-ray tube target	Cluster microfocal spot X-ray tube, Cu or Co, Cr, Fe, Ni, Mo, Ag, W target (optional according to test sample)
X-ray tube voltage	Maximum 50kV, 0-50kV adjustable
X ray tube power	Maximum 40W, 0-40W adjustable
XRF energy resolution	127eV@8keV
XRF detection range	Mg—U
Cooling Method	Water Cooling
Sample particle size	Sample particle <150um (100 mesh sieve)
Sample weight	about 20mg
Working temperature	-10°C-35°C
Weight	15KG
The power supply	Lithium ion battery or power adapter
Size	50×40×18.8cm□ L×W×H□
CrystalX analysis software language	English and other languages



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

