



DSC DIFFERENTIAL SCANNING CALORIMETER - TORONDSC™

The ToronDSC™ Differential Scanning Calorimeter is designed for fast and accurate measurement of heat flow and specific heat across a wide range of materials. It provides reliable thermal analysis for both routine quality control and advanced research applications. With a stable metal furnace and a robust system design, ToronDSC™ delivers consistent results that support confident material evaluation.

Built for efficiency and dependable operation, ToronDSC™ combines dual atmosphere capability with real-time monitoring to ensure stable testing conditions throughout each run. The system requires minimal calibration, helping reduce setup time while maintaining high measurement accuracy. Its balanced design makes complex thermal analysis easier to manage, even in busy laboratory environments.

ToronDSC™ supports precise evaluation of material behavior during thermal transitions. It is widely used in materials research, polymer analysis, formulation development, and performance validation. Designed for laboratories that value accuracy without unnecessary complexity, ToronDSC™ offers a practical solution for both daily testing and in-depth thermal studies.

Key Highlights

- Fast and accurate DSC heat flow analysis
- Stable metal furnace for consistent thermal performance
- Dual atmosphere operation for flexible testing conditions
- Real-time monitoring to improve workflow efficiency
- Minimal calibration requirements for reliable daily use
- Suitable for quality control labs and R&D facilities

APPLICATIONS

DSC Differential Scanning Calorimeter – ToronDSC

DSC Differential Scanning Calorimeter – ToronDSC

Applications of DSC Differential Scanning Calorimetry

DSC supports laboratories that need to measure thermal transitions, compare formulations, and confirm material stability under controlled heating or cooling conditions.

Polymers and Plastics

Measure glass transition, crystallization, melting behavior, and oxidation onset to compare grades, optimize processing windows, and evaluate thermal history.

Pharmaceuticals

Study polymorphism, purity, compatibility, and stability of active ingredients and formulations during development, quality control, and storage studies.

Adhesives, Resins, and Composites

Evaluate curing reactions, crosslink density, and batch-to-batch consistency for thermosets, prepregs, and formulated materials.

Research and Failure Analysis

Investigate unexpected thermal events, compare unknown materials, and correlate heat-flow behavior with performance issues in development and troubleshooting work.

Key Capabilities of ToronDSC

ToronDSC analyzes thermal properties and detects critical transitions such as glass transition temperature (T_g), cold crystallization, melting behavior, oxidation onset, curing reactions, and broader phase changes that affect product stability and processing.

Understanding the Differential Scanning Calorimetry Curve

The DSC curve shows how heat flow changes as the sample is heated or cooled relative to a reference. Peak direction, onset temperature, and peak area help identify endothermic and exothermic events and quantify the heat involved in each transition.

Differential Scanning Calorimetry Curve

- I. Secondary transition:** a small baseline shift often linked to structural rearrangement or a minor glass transition.
- II. Endothermic peak:** heat absorbed during events such as melting or fusion.
- III. Additional endothermic event:** heat absorbed during decomposition or bond-cleavage reactions.
- IV. Exothermic peak:** heat released during crystallization, curing, or crosslinking.

FEATURES

1. Advanced Metal Furnace Design

The newly engineered metal furnace structure provides a more stable baseline and improved measurement precision. Indirect conduction heating ensures superior uniformity and stability, significantly reducing pulse radiation compared to traditional heating methods.

2. Dual Atmosphere Flow System

Features automatic switching between two-way atmosphere flows with rapid transition and minimal stabilization time. Includes an additional protective gas inlet for enhanced safety and flexibility.

3. USB Communication Interface

Equipped with a versatile and reliable USB interface for seamless, uninterrupted communication. Supports automatic reconnection for added convenience.

4. Programmable Multi-Stage Temperature Control

Allows for fully automated temperature programming with precise control over each stage.

5. 7-Inch Industrial-Grade Color Touchscreen

High-resolution display provides clear, comprehensive data visualization and an intuitive user interface for easy operation.

6. Enhanced Sensitivity and Accuracy

Upgraded components significantly increase measurement sensitivity and accuracy. The bidirectional operating system enables real-time spectral acquisition and online data analysis through integrated software.

THEORY & METHOD

Differential Scanning Calorimetry (DSC) is a thermal analysis technique used to measure the difference in heat flow between a test sample and a reference as both are subjected to a controlled temperature program. This method provides valuable insight into material behavior during physical or chemical transitions, such as melting, crystallization, or curing.

The DSC instrument used here is a heat-flow-type DSC calorimeter, which records the thermal flow difference in milliwatts. These measurements reflect the enthalpy characteristics of the material—specifically, heat absorbed or released during transitions.

DSC Curve Interpretation

In the diagram above:

1. The **x-axis** represents temperature (T) or time (t), increasing from left to right.
2. The **curve labeled T** shows the reference temperature profile, while the **AH-line** curve displays the sample's deviation.

When the sample and reference have nearly equal heat capacities and no thermal events occur, the temperature difference (ΔT) between them remains near zero, resulting in a **flat baseline**, as seen in regions **AB, DE, and GH**.

Once a thermal event begins, the sample's behavior deviates:

1. **C-D segment** shows an **endothermic dip**, typically due to melting or decomposition, where the sample absorbs heat.
2. **E-F-G segment** features an **exothermic peak**, indicating heat release, such as crystallization or chemical crosslinking.

Each peak's position, height, width, and direction reveals critical information about the sample's thermal behavior. The area under the peak correlates with the magnitude of the heat involved in the transition. These markers help identify the type of thermal event (endothermic vs. exothermic), its onset temperature, and the transition range.

While ideal DSC curves are often smooth and clearly defined, real-world results may be more complex due to the kinetic factors and reaction rates of the material. The testing environment, sample preparation, and heating rate also influence the curve's accuracy and reproducibility.

This detailed curve interpretation is essential in industries such as polymer development, pharmaceutical testing, and materials research, where thermal properties directly impact product performance and stability.

TECHNICAL SPECIFICATIONS

Technical Parameters	ToronDSC™ -100A	ToronDSC™ -300	ToronDSC™ -300C
Temperature Range	Room temperature~600°C		-40~600°C
Temperature Resolution	0.01°C	0.001°C	0.001°C
Temperature Fluctuation	±0.1°C	±0.001°C	±0.001°C
Temperature Repeatability	±0.1°C	±0.01°C	
Heating Rate/Cooling Rate	0.1~100°C/min		0.1~100°C/min/0.1~40°C/min
Temperature Control Method	PID temperature control, temperature rise, Constant temperature	PID temperature control, heating, constant temperature, cooling	
DSC Range	0~±600mW		0~±800mW
DSC Accuracy	0.01mW		
Voltage	AC220V/50Hz or customized		
Gas Flow Rate	0~300mL/min		
Gas Pressure	≤5mPa		

Technical Parameters	ToronDSC™ -100A	ToronDSC™ -300	ToronDSC™ -300C
Display Method	24-bit color, 7-inch LCD touch screen display		
Data Interface	Standard USB interface		
Product Photos	DSC Differential Scanning Calorimeter – ToronDSC	DSC Differential Scanning Calorimeter – ToronDSC	DSC Differential Scanning Calorimeter – ToronDSC

Instrument Accessories

Variety of crucibles are available, including ceramic crucibles, aluminum crucibles, and crucibles with lids.

DSC Differential Scanning Calorimeter – ToronDSC

Tablet Press

Instrument Interface

The Initial Status key displays ambient temperature, sample temperature, and other relevant parameters.

DSC Differential Scanning Calorimeter – ToronDSC

The Parameter Settings key is used to configure experimental parameters, typically through the software interface.

DSC Differential Scanning Calorimeter – ToronDSC

The Device Information key displays detailed information about the instrument. The Administrator Channel is used by authorized personnel for internal temperature calibration.

DSC Differential Scanning Calorimeter – ToronDSC

The Start Run key displays real-time data after the operation is initiated through the computer software.

Test Graph

DSC Differential Scanning Calorimeter – ToronDSC

DSC Differential Scanning Calorimeter – ToronDSC



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

