



# **3 ZONES ROTARY TILTING TUBE FURNACE 1200°C -TT-RTF-120- 1200 PRO**

The 3 Zones Rotary Tilting Tube Furnace 1200°C - TT-RTF-120-1200 Pro is built for advanced thermal processing, integrating an automated rotation system, feeding and discharging functions, and a tilt-adjustable furnace body. Its furnace tube rotates a full 360° with adjustable speed to ensure uniform heating, while the feeding speed can be fine-tuned in real time for precise control. The body tilts up to 15° for easier material handling, and a collection unit at the tube's outlet enables continuous small-scale production, making it well-suited for research and pilot-scale testing.

This model comes with a  $\Phi 120 \times 1200$  quartz tube and three independent heating zones (200 mm each), offering stable continuous operation at 1100°C and a maximum temperature of 1200°C. Powered by 220V/13 kW (with a 110V option available), it utilizes high-quality alloy resistance wire (HRE) heating elements and an N-type thermocouple, ensuring accurate temperature control, efficiency, and durability for laboratory and industrial applications.

## APPLICATIONS

---

3 Zones Rotary Tilting Tube Furnace 1200°C -TT-RTF-120-1200 Pro

3 Zones Rotary Tilting Tube Furnace 1200°C -TT-RTF-120-1200 Pro

**Advanced Ceramic Materials** – Enables precise sintering, calcination, and densification of ceramics. Delivers uniform heating across three zones, ensuring superior mechanical strength and thermal stability for both research and industrial applications.

**Nanomaterials & Functional Composites** – Provides accurate, programmable thermal environments for synthesizing and modifying nanostructures, catalysts, and composites. Supports R&D in energy storage, electronics, and advanced coatings.

**Semiconductor Processing** – Facilitates high-temperature annealing, oxidation, and diffusion treatments for wafers and semiconductor substrates. Offers precise zone-controlled heating critical for microelectronics and optoelectronics.

**Powder Metallurgy** – Supports continuous small-scale sintering of metal powders into dense, uniform parts. The integrated feeding, rotation, and collection system makes it ideal for pilot production and alloy development.

**Petrochemical & Catalyst Research** – Suitable for catalyst activation, reforming studies, and thermal decomposition under controlled atmosphere or vacuum. Ensures reproducible performance for chemical laboratories.

**Aerospace Materials Engineering** – Enables processing of high-temperature alloys, ceramics, and composites used in aerospace applications. Validates thermal stability and mechanical strength under extreme thermal conditions.

**Energy & Battery Materials** – Provides thermal treatment for solid-state electrolytes, lithium-ion battery cathodes, and next-generation energy storage materials. Ensures consistency, purity, and high-temperature reliability.

**Universities & Research Institutes** – Ideal for multidisciplinary studies in physics, chemistry, materials science, and engineering. Combines automation, safety, and precision for advanced academic experiments.

**Pilot-Scale Industrial Applications** – Equipped with automatic feeding and discharging systems, it supports continuous small-batch processing for R&D-to-production transition. Perfect for scale-up testing and certification purposes.

## FEATURES

---

### Furnace Shell

- The shell is constructed from Q235 low-carbon steel with a corrosion-resistant, electrostatically coated surface. A durable double-layer metal frame is integrated with a heat-insulated fan structure to maintain a low external surface temperature. Designed for horizontal operation, the furnace features a top-opening configuration for easy access.

### Refractory Material

- The furnace includes a multi-layer insulation system made of lightweight alumina ceramic fiber and high-grade insulation board. It contains no asbestos and offers excellent thermal efficiency with minimal heat loss and reduced energy consumption.

Heating System

- Premium HRE spiral resistance wire heating elements are securely embedded in the furnace’s inner walls, delivering radiant heat throughout the chamber for energy-efficient and long-lasting performance.
- Ring-type resistance wire heating ensures consistent temperature distribution.

Temperature Control Panel

- Equipped with advanced controllers designed for complex process requirements.
- Supports up to 30 programmable segments for automatic temperature curve control.
- Features over-temperature and fault alarms, along with automatic safety shutoff.

Rotary and Automatic Feeding System

- The furnace tube uses a gear-driven electric motor to rotate 360° at an adjustable speed of 12-20 rpm, ensuring even heat exposure. The furnace body is electrically tiltable, allowing for smooth material entry and discharge. An automatic screw feeder with adjustable speed delivers consistent material flow, and a collection unit at the tube outlet enables continuous processing.

TECHNICAL SPECIFICATIONS

|              |                 |                  |                  |                  |
|--------------|-----------------|------------------|------------------|------------------|
| Item         | TT-RTF-60-12PRO | TT-RTF-100-12PRO | TT-RTF-150-12PRO | TT-RTF-200-12PRO |
| Tube Size    | Φ60×1200mm      | Φ100×1200mm      | Φ150×1500mm      | Φ200×1800mm      |
| Heating Zone | 440mm           | 440mm            | 600mm            | 800mm            |

|                        |                                                                                                                                                                                                                                          |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Max Temperature        | 1200°C                                                                                                                                                                                                                                   |
| Work Temperature       | 1100°C                                                                                                                                                                                                                                   |
| Heating Element        | HRE high temperature alloy resistance wire                                                                                                                                                                                               |
| Chamber Material       | Polycrystalline inorganic alumina ceramic fiber material prepared by wet vacuum filtration                                                                                                                                               |
| Temperature accuracy   | ±1°C                                                                                                                                                                                                                                     |
| Thermocouple           | N type thermocouple                                                                                                                                                                                                                      |
| Temperature Controller | Intelligent microcomputer PID temperature control instrument, SSR voltage regulation control, PID parameter self-tuning function; programmable 32 program segments, program temperature rise, program heat preservation, program cooling |
| Heating Rate           | 1-20°C/min freely adjustment                                                                                                                                                                                                             |
| Furnace Structure      | Furnace temperature control integrated structure                                                                                                                                                                                         |
| Protection             | Modular control, sound and light alarm signal will be sent out for over temperature and broken coupler in the working process, and the protection action will be completed automatically                                                 |
| Safety Protection      | The equipment is equipped with a circuit breaker, which will automatically pop open in case of short circuit leakage, which can protect the safety of equipment and operators                                                            |
| Furnace Shell          | High-quality cold-rolled steel plate CNC machine tool blanking processing, after welding, grinding, polishing, phosphating, pickling, surface electrostatic spraying plastic powder                                                      |
| Certificate            | ISO9001 CE SGS TUV                                                                                                                                                                                                                       |
| Standard Accessories   | 1□ Furnace Body 1 set 2□ Power cable 3 meter 3□ Thermocouple 1 pc 4□ Manual 1 pc 5□ Furnace Hook 1 pair<br>6□ High Temperature 1 pair                                                                                                    |
| OEM                    | According to the output requirements of the actual material of the user, the customization of special specifications of the product can be realized                                                                                      |



**Torontech Inc.**

251 Consumers Road, Suite 1200  
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652  
Email: [info@torontech.com](mailto:info@torontech.com)

**Torontech, LLC**

601 Brickell Key Drive, Suite 700  
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919  
Email: [sales@torontech.com](mailto:sales@torontech.com)

**Torontech ME FZE**

Meydan Grandstand, 6th Floor,  
Meydan Road, Nad Al Sheba,  
Dubai-261440, UAE 

Phone: +971 488 19252  
Email: [me@torontech.com](mailto:me@torontech.com)

