

Dual Drum Friability Tester- ToronPharma-F2

Features:

- Complies with USP, IP, Ph. Eur. specifications
- Patented automatic sample discharge
- Loading of samples into the drum without removing the drum
- Unique front loading system of the drum
- 10o tilt facility as per USP
- Supports standard and abrasion drum
- Programmable counter and timer
- Snap-click knob to hold the drum



Specifications

Speed	25 RPM fixed
Accuracy	+1 RPM
Timer	Settable from 1 sec. to 9 hrs. 59 mins. 59 secs.
Counter	1 to 9999 revolutions
No. of Drums	2
Power Supply	110 V.A.C., 60 Hz, 20 VA
Display	5 digit LED for elapsed time & rev. count
Dimension (L x W x H)	(350 x 310 x 430) mm approx.
Weight	Net 13 kgs- Gross 18 kgs

Optional accessory:

- AD drum™ (Qty. 2)
- Removable trays (Qty. 2)

IQ / OQ documents are available

Part #	Description
ToronPharma-F2	Dual Drum Friability Tester
	Features:
	<ul style="list-style-type: none"> • Complies with USP, IP, Ph. Eur. specifications • Patented automatic sample discharge • Loading of samples into the drum without removing the drum • Unique front loading system of the drum • 10o tilt facility as per USP • Supports standard and abrasion drum • Programmable counter and timer • Snap-click knob to hold the drum
	Standard Accessory Included:
	<ul style="list-style-type: none"> • AD drum™ (Qty. 2) • Removable trays (Qty. 2)
TT-006	IQ / OQ documents



<http://www.torontech.com/pharmaceutical/friability-testers>

Torontech Friability Testers



Dual Drum Friability Tester- ToronPharma-F2W

Features:

- Complies with USP, IP, Ph. Eur. specifications
- Patented automated sample discharge
- Variable RPM (20 to 50)
- Dual drum
- Programmable Counter and Timer
- Unique front loading system of the drum
- Results like % friability and weight-loss can be calculated
- Interfaced to Mettler, Sartorius and Contech weighing scales
- Printout of the results is possible
- Snap-click knob to hold the drum
- 10o tilt of the drum as per USP
- Supports standard or abrasion drum
- Online validation of RPM



Specifications

Speed	Variable speed from 20 to 50 RPM
Accuracy	+0.5 RPM
Timer	Settable from 1 sec. to 9 hrs. 59 mins. 59 secs.
Counter	1 to 9999 revolutions
No. of Drums	1
Power Supply	100/AC, 50/60 Hz, 550 VA
Display	24 x 2 character LCD
Dimension (L x W x H)	(350 x 260 x 430) mm (approx.)
Weight	Net 11 kgs- Gross 15 kgs

Standard Accessory Included:

- AD drum™ (Qty. 2)
- Removable trays (Qty. 2)

IQ / OQ documents are available

Part #	Description
ToronPharma-F2W	Dual Drum Friability Tester
	Features:
	<ul style="list-style-type: none"> • Complies with USP, IP, Ph. Eur. specifications • Patented automated sample discharge • Variable RPM (20 to 50) • Dual drum • Programmable Counter and Timer • Unique front loading system of the drum • Results like % friability and weight-loss can be calculated • Interfaced to Mettler, Sartorius and Contech weighing scales • Printout of the results is possible • Snap-click knob to hold the drum • 10o tilt of the drum as per USP • Supports standard or abrasion drum • Online validation of RPM
	Standard Accessory Included:
	<ul style="list-style-type: none"> • AD drum™ (Qty. 2) • Removable trays (Qty. 2)
TT-007	IQ / OQ documents



<http://www.torontech.com/pharmaceutical/friability-testers>