

## Manual Powder Flow Tester – ToronPharma-FT01

### Features:

- Complies with Ph. Eur. Specifications
- Provides simple, repeatable flowability index of
- Powders and granules
- Easy calculation of angle of repose
- Outflow openings- 10 mm, 15 mm and 25 mm as per Ph. Eur.
- Customized openings – 5 mm, 10 mm, 15 mm etc. are also available



Part #	Description
ToronPharma-FT01	Manual Powder Flow Tester
	<b>Features:</b>
	<ul style="list-style-type: none"> <li>• Complies with Ph. Eur. Specifications</li> <li>• Provides simple, repeatable flowability index of</li> <li>• Powders and granules</li> <li>• Easy calculation of angle of repose</li> <li>• Outflow openings- 10 mm, 15 mm and 25 mm as per Ph. Eur.</li> <li>• Customized openings – 5 mm, 10 mm, 15 mm etc. are also available</li> </ul>
TT-004	<b>IQ/OQ documents</b>

<http://www.torontech.com/pharmaceutical/manual-powder-flow-tester-eft>