



ICU BEDS

Electric ICU beds built to IEC 60601-2-52, from multifunction to column-lift.

ICU beds give intensive care teams full electric control of patient positioning, from backrest and leg articulation to Trendelenburg tilt and full-height travel. The ToronCare™ ICU range from Torontech is built to IEC 60601-2-52, the international standard for medical beds, and serves hospitals, surgical recovery units, and high-dependency wards.

Every model in the range carries a 240 kg safe working load, an IP54 ingress rating, and electric and mechanical CPR for emergency response. The four configurations differ in weighing, lateral tilt, and lift design, so you can match the bed to the acuity and workflow of each unit.

Applications

ToronCare™ ICU beds support acute and long-term care settings where positioning, monitoring, and rapid intervention shape patient outcomes.


- Intensive Care Units (ICU): continuous positioning, Trendelenburg tilt, and CPR access for ventilated and immobile patients.
- Cardiac Care Units (CCU): one-button cardiac chair positioning to seat patients upright and ease breathing.
- Dialysis and Renal Wards: in-bed weighing on scale-equipped models to track fluid removal without patient transfers.
- Pressure Injury and Wound Care: lateral tilt to offload bony prominences and assist turning.
- Emergency and Post-Anesthesia Units: rapid flat positioning and height adjustment for resuscitation and recovery.
- Long-Term Acute Care: durable surfaces, central braking, and low minimum heights for extended stays and fall prevention.

Why Choose Torontech for Your ICU Bed

Torontech North America supplies testing and measurement equipment to laboratories, quality teams, and manufacturers, backed by responsive technical support and calibration service. The ICU bed range is built for reliable, standards-driven results, so you can source confidently from a single partner. Our team helps you match the right configuration to your samples, applicable standards, and throughput, then supports you through installation and beyond.

Request a Quote

Tell us the acuity of your unit and whether you need weighing, lateral tilt, or a column lift, and we will recommend the right ToronCare™ ICU bed. [Contact the Torontech](#) team for pricing and lead times.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

