

HEMOCARE HOSPITAL BED - HEMOCARE BED



Homecare Hospital Beds | Manual, Electric & Hi-Low Models to IEC 60601-2-52

A homecare bed brings hospital-grade positioning into the home, helping caregivers reposition, lift, and care for bedridden or recovering patients safely. The ToronCare line covers manual and electric homecare hospital beds built to IEC 60601-2-52, from simple fixed-height frames to fully automated nursing beds. Each model suits a different level of care, mobility, and budget.

Families, nursing homes, and long-term care facilities use these beds to lower caregiver strain, prevent pressure ulcers, and keep patients comfortable for long periods. This range spans manual cranks, push-button electric controls, hi-low height adjustment, lateral turning, and built-in toileting, so you can match the bed to the patient and the setting.

Applications

Homecare hospital beds support extended care wherever a patient spends long periods in bed, at home or in a facility.

- Home-Based Long-Term Care: gives families a safe, adjustable bed for relatives with chronic illness or limited mobility.
- Post-Surgical and Rehabilitation Recovery: backrest and leg positioning help patients rest, sit up, and recover after surgery or injury.
- Elderly and Dementia Care: low heights, fall-prevention rails, and stable bases protect aging and at-risk patients.
- Stroke and Neurological Recovery: turning and tilt functions reposition patients who cannot move themselves.
- Palliative and Hospice Care: gentle positioning and in-bed hygiene keep patients comfortable and dignified.
- Nursing Homes and Long-Term Care Facilities: durable frames and simple controls support continuous, multi-shift use.

Homecare Hospital Bed Range

The ToronCare range runs from entry-level fixed-height beds to a premium automated model. Use the table to compare how each bed operates and who it suits.

Model	Type & Operation	Description
ToronCare™ 1080	Electric, stationary (no castors)	Stationary bed on a four-piece stand-board base instead of castors, with remote backrest, leg-rest, and hi-low height. The non-rolling design suits fall-risk and cognitively impaired patients who need a stable, stay-put bed.
ToronCare™ 5081	Manual, multifunction	Hand-cranked bed that needs no electricity, with lateral turning, a chair position, a built-in toilet opening, and a bed-side shampoo station. A budget choice for off-grid homes and cost-conscious long-term care.
ToronCare™ 5083	Electric, multifunction	Push-button bed with electric lateral turning, a one-key chair position, electric toileting, a one-button CPR flat, and a back-up battery for outages. Suits caregiver-intensive and emergency-ready home care.
ToronCare™ 1085	Electric, foldable hi-low	Cost-effective folding bed that one person can assemble without tools, with hi-low height (335 to 725 mm), backrest, knee bend, and Trendelenburg tilt. Easy to ship, store, and move between rooms.
ToronCare™ 5086	Electric, hi-low adjustable	Heavy-duty bed with a wide 240 to 640 mm height range, whole-bed Trendelenburg tilt, and a 240 kg capacity. The ultra-low height protects fall-risk patients, and an optional stand-assist rail aids transfers.
ToronCare™ 5088	Electric, premium automated	Flagship bed with automated timed turning, voice-guided controls, an anti-slide chair, a powered commode pan, and an emergency stop. Built for high-dependency and around-the-clock care.
ToronCare™ 1090	Electric, fixed-height	The simplest and lightest model at 80 kg, with a fixed 380 mm height and remote backrest (to 75°) and leg-rest control. A straightforward, low-maintenance choice for general and outpatient care.

How to Choose the Right Homecare Bed

- For the simplest, lowest-maintenance setup: choose the fixed-height ToronCare™ 1090, which covers backrest and leg positioning without height or turning mechanisms.
- For a stable, non-rolling bed for fall-risk or dementia patients: choose the stationary ToronCare™ 1080 on its stand-board base.
- For easy shipping, storage, and one-person setup: choose the foldable ToronCare™ 1085, which adds hi-low height and Trendelenburg tilt.
- For care without reliable electricity or on a tight budget: choose the manual ToronCare™ 5081, with turning, chair, toilet, and shampoo functions on hand cranks.

- For powered multifunction care with emergency readiness: choose the electric ToronCare™ 5083, with electric turning, a CPR flat, and a back-up battery.
- For heavy-duty, ultra-low fall protection with tilt: choose the ToronCare™ 5086, with a 240 kg capacity and a 240 to 640 mm height range.
- For high-dependency, around-the-clock care: choose the premium ToronCare™ 5088, with automated timed turning, voice prompts, and a powered commode.

Why Choose Torontech for Your Homecare Bed

Torontech supplies testing and measurement equipment to laboratories and manufacturers across North America and worldwide, backed by responsive technical support and calibration service. The homecare bed range is built for reliable, standards-driven results, so you can standardize on one supplier as your testing needs grow. Our team helps you match the right configuration to your samples, applicable standards, and throughput, then supports you through installation and beyond.

Request a Quote

Tell us the patient's mobility, the level of care, and your budget, and we will point you to the right ToronCare homecare bed. [Contact the Torontech team](#) for pricing, lead times, and the nearest distributor.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

