



HOSPITAL ICU/CCU BED - TORONCARE™ 1076

Hospital ICU/CCU Beds - ToronCare™ 1076 is our new designed hospital bed to meet multi-function requirements of patients.


APPLICATIONS

- **Intensive and Critical Care Units** – Provides multifunctional support with adjustable back, knee, and Trendelenburg positions.
- **Emergency and Post-Surgical Care** – Features quick CPR positioning, battery backup, and weighing scale for continuous monitoring.
- **Cardiac and Respiratory Support** – Chair function assists with lung ventilation and patient comfort.
- **Long-Term Patient Care** – Reduces fall risk with protective side rails, exit alarms, and angle indicators.
- **Diagnostic and Imaging Use** – X-ray translucent backboard allows imaging without moving the patient.
- **Hospital Infection Control** – Easy-to-clean electro-coated and powder-coated steel design ensures hygienic use.
- **Safe Mobility and Handling** – Equipped with silent wheels, central brake, and bed extension for patient stability and convenience.



TECHNICAL SPECIFICATIONS

Technical Parameters	Other Features		
Outer Dimensions	2200 x 1050 mm	Chair Function	Battery backup
Height	525 - 965 mm	CPR Function	Weighing scale
Length Extension	200 mm	Footboard curve function	4 double-side silent wheels
Back Section	≥ 75°	Treated by electro-coating and powder coating	All screws made of stainless steel material
Knee Section	≥45°	Fall Protection	Angle indicator panel
Trend./Reverse Trend.	≥15°	Cold-rolled steel bed boards	Height adjustable 304 stainless
Scale	0.1 kg precision	PP material handrails	Drainage hook
Off bed Alarm	Together with scale	Side Rail controls	With anti-bumper, lighting lamp
Working Load	300 kg	Bed Extension	Central-controlled side brake
Net Weight	160 kg	X-ray backboard	6pc Linak Motors
Gross Weight	200 kg	Head and Foot Boards removable with Locks	One step molding powder coated steel dust cover
Package size	2200 x 900 x 60 mm	Embedded operator	Content



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

