



PREMIUM ELECTRIC HOSPITAL BED - TORONCARE™ 1060

Critical care demands rapid response capabilities. The ToronCare™ 1060 answers with a quick-release CPR lever that instantly flattens the back-rest during emergencies, complemented by backup battery protection and reset positioning. This premium electric hospital bed combines touchpad controls with custom bed board engineering that enhances patient comfort through improved flexibility. At 148 kg with 2260 x 1050 mm dimensions, it supports up to 250 kg while delivering professional-grade reliability for intensive care units and post-surgical recovery.

Simultaneous back-rest and leg-rest adjustment streamlines caregiver workflow, while Trendelenburg positioning reaches 10° for therapeutic interventions. Electric central brake castors eliminate manual effort, and gas spring-activated drop-down railings provide secure patient protection. Touchpad interface with indicator lights replaces failure-prone mechanical keypads for long-term durability. Height adjusts from 530-710 mm, back-rest articulates to 65°, and leg-rest positions to 40°. Optional X-ray radiolucent sections enable bedside imaging, while I.V pole brackets at all four corners maximize treatment flexibility.

Rehabilitation facilities and general hospital wards benefit from this adjustable medical bed's robust steel and ABS & PP construction, with detachable head and foot boards simplifying maintenance in demanding healthcare environments.

APPLICATIONS

- **Critical and Intensive Care** – Offers electrically controlled adjustments, emergency CPR function, and reset positioning for quick response.
- **Post-Surgical Recovery** – Provides smooth adjustment of back-rest and leg-rest simultaneously for patient comfort.
- **General Hospital Wards** – Ensures reliable care with durable construction, protective side rails, and electric central braking.
- **Emergency Preparedness** – Backup battery and quick-release CPR lever allow safe operation during power outages or urgent situations.
- **Diagnostic Support** – Optional X-ray radiolucent sections enable imaging without moving the patient.
- **Long-Term and Rehabilitation Care** – Enhances safety, comfort, and mobility with adjustable functions and ergonomic design.

FEATURES

- Built with durable steel framework and high-quality ABS & PP plastic
- Controlled by user-friendly touchpad (replacing traditional keypads, which have limited service life due to mechanical failure)
- Custom bed board improves softness and bending of the bed, increasing the comfort of patients
- Electrically controlled bed positions provide easy adjustment of patient posture, with an indicator light to show the work status
- Back-rest and leg-rest can be adjusted up/down at the same time for more efficient operation
- Quick release lever on back-rest for emergency CPR
- Reset button adjusts the bed to the flat / lowest setting from any position
- Backup battery provides protection for emergency power failure
- Detachable head and foot board made from high-quality PP plastic
- 4 pcs drop-down PP protective railing, controlled by gas spring
- Central brake castors controlled by electric operation
- IV pole mount bracket at all 4 corners of the bed frame

Standard Accessories (Included)

- IV Pole

Optional Accessories (Not Included)

- Monkey pole
- X-ray radiolucent head and back section

TECHNICAL SPECIFICATIONS

Dimensions	
Overall dimensions (L x W x H)	2260 mm x 1050 mm x 530-710 mm (89.0 in x 41.3 in x 20.9-28.0 in)
Hi-Low	710 mm-530 mm (28.0-20.9 in)
Bed Board dimensions (L x W)	1940 mm x 900 mm (76.4 in x 35.4 in)
Diameter of universal central locking castors	150 mm (5.9in)
Bed Articulations	
Back-rest	0-65°
Leg-rest	0-40°
Trendelenburg and Rev. Trendelenburg	0-10°
Weight	
Net Weight	148 Kg (326 lbs)
Load Capacity	<250 Kg (<551 lbs)



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

