



ELECTRIC HOSPITAL BED - TORONCARE™ 1040

The ToronCare™ 1040 is a professional three-function electric hospital bed engineered for hospital wards, post-surgical recovery, rehabilitation support, and long-term care facilities. Constructed with robust steel framework and premium ABS & PP plastic, this adjustable medical bed provides reliable patient positioning with a 250 kg load capacity and professional-grade functionality for home healthcare.

This electric adjustable hospital bed features remote-controlled back-rest (0-70°), leg-rest (0-35°), and hi-low adjustment (445-655 mm) for enhanced patient comfort and therapeutic positioning. The conveniently located brake pedal system with 150 mm central castors ensures secure mobility and stability during patient care. With dimensions of 2160 x 1150 mm and weighing 130 kg, it includes a standard I.V pole for immediate clinical use.

Ideal for extended patient stays and rehabilitation therapy, the ToronCare™ 1040 improves caregiver accessibility through ergonomic positioning and intuitive controls. This durable medical bed delivers dependable performance for everyday clinical environments and home-based care applications.

APPLICATIONS

- **General Patient Care** – Provides essential electric adjustments for back-rest, leg-rest, and bed height.
- **Post-Surgical Recovery** – Facilitates safe and comfortable positioning during the healing process.
- **Long-Term Care Facilities** – Enhances daily comfort and caregiver accessibility for extended patient stays.
- **Rehabilitation Support** – Assists with mobility and therapy needs through ergonomic positioning.
- **Home Healthcare** – Offers professional-grade functionality for patients requiring care at home.
- **Hospital Wards** – Durable design ensures reliable use in everyday clinical environments.

FEATURES

- Built with durable steel framework and high-quality ABS & PP plastic

- Three functions easily adjusted by remote control (back-rest, leg-rest and high-low)
- Brake pedal is conveniently located below the foot board

Standard Accessories (Included)

- I.V Pole

TECHNICAL SPECIFICATIONS

Dimensions	
Overall dimensions (L x W x H)	2160 mm x 1150 mm x 445-655 mm (85.0 in x 45.3 in x 17.5-25.8 in)
Hi-Low	655 mm - 445 mm (25.8 - 17.5 in)
Bed Board dimensions (L x W)	1935 mm x 900 mm (76.2 in x 35.4 in)
Diameter of central brake castors	150 mm (5.9 in)
Bed Articulations	
Back-rest	0-70°
Leg-rest	0-35°
Weight	
Net Weight	130 Kg (287 lbs)
Load Capacity	<250 Kg (<551 lbs)



TORONTECH

Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com



USA | Canada | EU | UAE | GCC | KSA | India | Africa | Latin America | Asia

