



ELECTRIC HOSPITAL BED - TORONCARE™ 1030

The ToronCare™ 1030 is a professional three-function electric hospital bed designed for hospital wards, rehabilitation facilities, post-surgical recovery, and home healthcare. Built with durable steel framework and quality ABS & PP plastic, this adjustable medical bed delivers reliable patient positioning with a 250 kg load capacity.

This electric adjustable hospital bed offers remote-controlled back-rest (0-70°), leg-rest (0-35°), and hi-low adjustment (445-655 mm) for optimal patient comfort and therapeutic positioning. The convenient brake pedal system with 150 mm central castors ensures secure movement and stability. With dimensions of 2160 x 1150 mm and including a standard I.V pole, it provides professional-grade functionality for clinical environments.

Ideal for long-term patient care and home-based medical settings, the ToronCare™ 1030 enhances caregiver ergonomics through intuitive controls and versatile positioning. This medical bed combines durability with user-friendly operation for everyday healthcare applications.

APPLICATIONS

- **General Hospital Wards** – Provides essential electric adjustments for improved patient care and comfort.
- **Post-Surgical Recovery** – Supports proper positioning of back-rest and leg-rest to aid healing.
- **Rehabilitation Facilities** – Assists patients with mobility and therapy through adjustable height and posture functions.
- **Long-Term Patient Care** – Ensures safety and ease of use for patients requiring extended stays.
- **Home Healthcare** – Delivers professional-grade electric functionality for patients under home-based care.
- **Clinical Use** – Durable and user-friendly design makes it suitable for everyday medical environments.

FEATURES

- Built with durable steel framework and high-quality ABS & PP plastic
- Three functions easily adjusted by remote control (back-rest, leg-rest and hi-low)

- Brake pedal is conveniently located below the foot board

Standard Accessories (Included)

- I.V Pole

TECHNICAL SPECIFICATIONS

Dimensions	
Overall dimensions (L x W x H)	2160 mm x 1150 mm x 445-655 mm (85 in x 45.3 in x 17.5-25.8 in)
Hi-Low	655 mm - 445 mm (25.8 - 17.5 in)
Bed Board dimensions (L x W)	1935 mm x 900 mm (76.2 in x 35.4 in)
Diameter of central brake castors	150 mm (5.9 in)
Bed Articulations	
Back-rest	0-70°
Leg-rest	0-35°
Weight	
Net Weight	112 Kg (247 lbs)
Load Capacity	<250 Kg (<551 lbs)



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

