

TORONDENT™ DENTAL UNITS



ToronDent™ 1000 series units are designed to make dental treatment sessions more comfortable and relaxing for your patients - making your work easier. With accurate positioning and support, ToronDent™ series units help you focus on providing high-quality clinical treatment. Torontech works closely with dental professionals and introduces smart new ways to ensure maximum comfort and exceptional quality treatment.

Dental Units Accessories

1. ToronDent™ 1100

Main Features

- Halogen LED Lamp
- Micro-controlled noiseless DC actuators
- Luxury cushions
- Position Memory function and assistance control
- Turnable Cuspidor
- Built in water purifying system
- Multifunction Foot Pedal
- Three way syringes

Technical Specifications

- Input Power: 220Vac, 50Hz or 110Vac, 60Hz
- Air Pressure: 0.5-0.8 MPa
- Water pressure: 0.2-0.4 MPa
- Weight: around 230kg (500 lbs.)

Standard Accessories

- 2 pcs. of high speed handpiece connector
- Low speed handpiece connector
- 2 pcs. of 3-way syringe (cold and hot)
- Operation round lamp
- X-ray Film viewer
- Saliva ejector (strong and weak)
- Automatic cup filler
- Rotatable ceramic spittoon
- Rotatable assistant device
- Built-in water purified system
- 24V D.C noiseless motor
- Sanitary cushion and backrest
- Separated air and water system
- Dentist Stool

Optional Accessories

- Hand Piece
 - Intraoral Camera
 - Ultrasonic Scaler
 - Curing Light
-

2. ToronDent™ 1200

Main Features

- Equipped with sterilization system of hoses
- Infrared sensor switch (non-touched), universal rotatable operating light
- Operating light handle removal option for sterilization

- Micro-controlled noiseless DC actuators
- Luxury cushions
- Position Memory function and assistance control
- Turnable Cuspidor
- Built in water purifying system
- Multifunction Foot Pedal

Technical Specifications

- Input Power: 220Vac, 50Hz or 110Vac, 60Hz
- Air Pressure: 0.5-0.8 MPa
- Water pressure: 0.2-0.4 MPa
- Weight: 230kg (500 lbs.)

Standard Accessories

- 2 pcs. of high speed handpiece connector
- Low speed handpiece connector
- 2 pcs. of 3-way syringe (cold and hot)
- Operation round lamp
- X-ray Film viewer
- Saliva ejector (strong and weak)
- Automatic cup filler
- Rotatable ceramic spittoon
- Rotatable assistant device
- Built-in water purified system
- 24V D.C noiseless motor
- Sanitary cushion and backrest
- Separated air and water system
- Dentist Stool

Optional Accessories

- Hand Piece
 - Intraoral Camera
 - Ultrasonic Scaler
 - Curing Light
-

3. ToronDent™ 1300

Main Features

- Separated air and water system
- High quality micro fiber cushion
- Various color variety options for cushion
- Luxury arm rest
- LED sensor light
- 90o rotatable cuspidor
- Multi-function foot pedal

Technical Specifications

- Input Power: 220Vac, 50Hz or 110Vac, 60Hz
- Air Pressure: 0.5MPa – 0.8 MPa
- Water Pressure: 0.2MPa-0.4MPa
- Weight: 240 Kg (529 lbs.)


Standard Accessories

- Computer controlled noiseless D.C motor

- Luxury micro fiber cushion and backrest
- LED sensor lamp (LED and halogen color)
- 3 memories function and assistant control system
- Turnable sandblasting glass cuspidor
- Built-in water purifying system
- Multifunction foot pedal
- Magnetic valve
- X-ray film viewer
- Strong and low suction
- 2 pcs. of 3-way syringe (cold and hot)
- Doctor stool

Optional Accessories

- Handpiece
- Intraoral camera set
- Ultrasonic scalar
- Curing light



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

