

LEAK TESTER LT-101P



With programmable vacuum levels, slow venting for improved results, validation report printing, and an oil-less vacuum pump, the Leak Tester LT-101P ensures reliable and efficient performance.

APPLICATIONS

- **Pharmaceutical Packaging** – Tests blister packs, vials, and bottles to ensure airtight sealing and product integrity.
- **Food and Beverage Industry** – Verifies leak-proof packaging to maintain freshness and prevent contamination.
- **Medical Device Manufacturing** – Confirms the integrity of sterile packaging for instruments and consumables.
- **Research and Quality Control Laboratories** – Provides reliable leak detection with programmable vacuum settings and validation reporting.
- **Industrial Packaging** – Ensures compliance and durability of sealed pouches, containers, and flexible packaging.
- **Regulatory Compliance Testing** – Supports validation in accordance with GMP and other international quality standards.

FEATURES

- Programmable vacuum level and hold time
- Slow venting for improved results
- Capability to print out validation reports
- Equipped with an oil-less vacuum pump for reliable performance.

TECHNICAL SPECIFICATIONS

Vacuum level	Better than 600 mm of Hg
--------------	--------------------------

Timer	Maximum can be set up to 99 minutes
	Minimum can be set up to 1 min
Vacuum setting level	Vacuum can stabilize to the set value +5 mm of Hg within 30 secs
Power requirement	230V. 50 Hz (110V is also available. Please specify when ordering)
Power	300 VA
Dimensions (L x W x H)	400 x 240 x 220 mm
Weight	Net: 17 kgs
	Gross: 21 kgs



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

