



HANDHELD XRF ANALYZER TRUEX SERIES

Light weight, detection range from Mg-U, with fast analysis speed, high accuracy, handheld and easy to carry, is the ideal tool for field testing.

The Handheld XRF Analyzer TrueX Series, a powerful handheld device designed for rapid and accurate analysis in various fields. With features like a compact design, advanced technology, and customizable reports, this Handheld XRF Analyzer ensures laboratory comparable results in the field.

It offers capabilities such as detecting key elements with high content and providing precise analysis without damaging the tested object. Additionally, its advanced Peltier electric refrigeration method eliminates the need for liquid nitrogen, reducing equipment complexity and maintenance costs. The TrueX series complies with national radiation safety standards, ensuring user safety and peace of mind.

1.5KG


Lightweight and easy to carry

25Elements

Up to 25 elements can be detected

DoubleBeam™

Built-in interlocking device, effectively protect user safety



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

