

SCOPEX CS DESKTOP SPECTROMETER FAMILY



Introducing the ScopeX CS Desktop Spectrometer Family, a cutting-edge range of spectrometers designed for high-precision analysis. These spectrometers feature advanced technology, including low-energy X-ray and intelligent vacuum systems, ensuring accurate detection of various elements. With options like Peltier electric refrigeration and automatic collimator and filter switching, the devices offer enhanced performance while reducing complexity and maintenance costs.

The fully enclosed metal chassis and leakage protection switch ensure radiation safety in compliance with international standards. Experience seamless data transmission through USB, WiFi, and Bluetooth connectivity, and customize test reports effortlessly with options for different file formats and personalization. With ScopeX CS Desktop Spectrometer Family, achieve precise analysis without compromising on quality or safety.

Reduce interference


Reduce spectral interference, thereby improving detection accuracy

Improve sensitivity

More accurate measurement of light elements

Wide range of applications

Suitable for steel, metallurgy and other fields



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

