



INSIGHT FAMILY (BENCHTOP)

The INSIGHT Family (Benchtop) – a versatile range of instruments designed to meet the diverse needs of industries such as electronics, semiconductor, and materials analysis. With its user-friendly operation and advanced features, INSIGHT ensures fast, accurate, and non-destructive analysis of coating thickness and composition, even for complex structures like PCBs and semiconductor wafers.

Equipped with intuitive software and automation capabilities, INSIGHT accelerates sample measurements, enhancing productivity and efficiency. From wafer plating to lead frame treatments, the INSIGHT Family (Benchtop) offers reliable solutions for your analytical needs.

Multi-collimator

Multiple collimators are optional and combinable, and the configuration is flexible

Nanoscale

Measurable precision down to nanometer level

Helium system

Optional helium system for excellent light element measurement



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

