



GOLD PURITY DENSITY TESTER

Gold Purity Density Tester - Density Meter | Densimeter | Specific Gravity Tester

Gold Purity Density Testers are equipped with advanced features tailored for precise measurement of gold purity and density.

Gold Purity Density Testers features an RS-232 & USB interface for seamless connectivity and infrared temperature sensors for automatic water temperature compensation. With quick and accurate determination of precious metal purity and gold Karat number in just two steps, it ensures reliable results without contamination or the need for sulfuric acid solutions.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

