



TRANSFORMER OIL TESTING EQUIPMENT

Transformer Oil Testing Equipment from Torontech supports insulating oil analysis, petroleum product testing, transformer maintenance, and electrical equipment quality control. This category includes instruments for flash point, moisture content, kinematic viscosity, oxidation stability, pour and cloud point, sulfur content, corrosive sulfur, density, dielectric breakdown voltage, petroleum color, and dissolved gas analysis.

Transformer oil testing helps power utilities, transformer service providers, petrochemical laboratories, and industrial maintenance teams evaluate oil condition, insulation performance, contamination risk, moisture level, oxidation behavior, gas formation, and low-temperature flow properties. These results support routine condition monitoring, oil acceptance, preventive maintenance, fault investigation, and long-term transformer reliability programs.

Torontech offers manual, semi-automatic, and automatic testing equipment for different laboratory workflows. Users can compare products by test parameter, sample type, automation level, diagnostic purpose, and testing requirement to select the right instrument for insulating oil analysis and petroleum quality control.

Compare Transformer Oil Testing Equipment by Testing Need

Use this table to match transformer oil and petroleum testing requirements with the right Torontech equipment.

Testing Need	Recommended Equipment	Purpose
Open cup flash point testing	Manual Cleveland Open Cup Flash Point Tester	Measures flash point behavior for transformer oil, lubricating oil, and petroleum product samples.
Automated open cup flash point testing	Automatic Cleveland Open Cup Flash Point Tester	Supports automated flash point testing for laboratories that need consistent operation and higher workflow efficiency.

Testing Need	Recommended Equipment	Purpose
Closed cup flash point testing	Manual Pensky-Martens Closed-Cup Flash Point Tester	Measures flash point using a closed-cup method for petroleum products and oil condition checks.
Automated closed cup flash point testing	Automatic Pensky-Martens Closed Cup Flash Point Tester	Provides automated closed-cup flash point testing for routine petroleum and transformer oil workflows.
Low-level moisture testing	Coulometric Karl-Fischer Titration Tester	Measures low water content in transformer oil, insulating liquids, and petroleum samples.
Moisture content testing	Volumetric Karl Fischer Moisture Titrator	Measures water content in oil and petroleum samples using volumetric moisture titration.
Kinematic viscosity testing	Manual Kinematic Viscosity Bath	Measures oil flow behavior under controlled temperature conditions for petroleum and insulating oil analysis.
Semi-automatic viscosity testing	Semi Automatic Kinematic Viscosity Apparatus	Supports routine kinematic viscosity testing with semi-automatic operation for lab workflows.
Low-temperature viscosity testing	Low Temperature Kinematic Viscosity Tester	Evaluates viscosity behavior at low temperatures to support cold-condition oil performance checks.
Oxidation stability testing	Oxidation Stability Tester for Mineral Insulating Oil	Evaluates oil aging behavior, oxidation resistance, and degradation tendency in mineral insulating oil.
Pour point and cloud point testing	Fully Automatic Pour and Cloud Point Tester	Measures low-temperature flow behavior and clouding characteristics of oil and petroleum samples.

Testing Need	Recommended Equipment	Purpose
Sulfur content analysis	X-ray Fluorescence Sulfur Content Apparatus	Measures sulfur content in petroleum products and related oil samples using XRF analysis.
Corrosive sulfur testing	Transformer Oil Corrosive Sulfur Tester	Evaluates corrosive sulfur tendency that may affect transformer oil condition and electrical components.
Density testing	Density Test by API Hydrometer	Measures density, relative density, or API gravity for petroleum and oil quality checks.
Dielectric breakdown testing	Insulating Oil Dielectric Breakdown Tester	Evaluates insulating oil strength under electrical stress and supports transformer oil condition monitoring.
Petroleum color testing	Colorimeter for Petroleum Products	Measures oil color for product quality review, condition checks, and petroleum sample comparison.
Dissolved gas analysis	GC for Analysis of Dissolved Gases in Transformer Oil (TOGA)	Supports transformer oil diagnostics by analyzing dissolved gases linked to electrical or thermal stress.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

