




MICROSCOPIC HYPERSPECTRAL IMAGING SYSTEM

Combining the benefits of a microscope and an imaging spectrometer, hyperspectral data acquisition of microscopic images can be performed at any time.

Microscopic Hyperspectral Imaging System can adapt existing biological microscopes, fluorescence microscopes, stereo microscopes, metallographic microscopes, and more, allowing them to be easily converted into hyperspectral microscopes. Users can customize microscope models based on their specific needs. The FigSpecR series of imaging spectrometers integrates both a visual camera and a hyperspectral camera. The visual camera allows for a quick preview of the sampled images, and once the images are confirmed to meet the requirements, hyperspectral image data collection can begin.

APPLICATIONS

System composition

Hyperspectral imaging spectroscopic camera (optional models: FS-20/FS-22/FS-23) x1, Lens x1, Microscope (any manufacturer's model can be specified) x1, PC application software x1

Applications

Microscopic Hyperspectral Imaging SystemMicroscopic Hyperspectral Imaging System



TORONTECH

Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com



USA | Canada | EU | UAE | GCC | KSA | India | Africa | Latin America | Asia

