



ELECTRIC HOSPITAL BED - TORONCARE™ 1085

APPLICATIONS

- **Homecare Settings** – Provides safe, adjustable support for patients receiving care at home with easy assembly and portability.
- **Long-Term Care Facilities** – Ensures comfort and reliability for daily patient needs with sturdy construction and safety rail protection.
- **Hospitals and Clinics** – Offers cost-effective electric adjustment functions suitable for general patient care.
- **Rehabilitation Centers** – Supports recovery by enabling customized positioning such as Trendelenburg and reverse Trendelenburg.
- **Elderly Care** – Enhances patient safety with 32 cm fall-protection side rails and ease of use for caregivers.
- **Emergency Preparedness** – Foldable design allows quick installation, efficient storage, and easy transport when needed.

FEATURES

- **Quality and Durability:** The homecare bed is sturdily constructed to match with the rigours of daily care. High-quality components make this a particularly durable, sturdy and low maintenance care bed that reduces further costs.
- **Cost-Effective:** This bed offers exceptional value for money, but not at the expense of materials or design.
- **Easy for Assembling & Carrying:** The homecare bed can be assembled, disassembled, and carried by one person without tools in few minutes.
- **Safety Side Fall Protection:** Mattresses or anti-decubitus mattress up to a thickness of 8 cm. The 32 cm siderail fall protection height is a great assurance for safety of users.

TECHNICAL SPECIFICATIONS

Functions	
Backrest Angle:	0-72°
Knee Section Angle:	0-40°
Height Adjust:	335 -725 mm
Trendelenburg and Reverse Trendelenburg angle	0~15°
Configurations	
Size:	2180 x 1000 x (335 -725) mm
Bed Board:	Wooden bar with high loading weight
Head and Foot Board:	Wooden
Guardrail:	Aluminum with surface processing, which is durable and warping preventing.
Castor	Universal castors with brakes
Optional Configuration	
Grip pole, Movable table, Bedside cabinet, Mattress	
Warranty	
Motor: 2 Years, Other Parts:1 year	
Specialty	
The bed is foldable which saves space for shipping and storage. It is easy to install the beds by one person without the help of any tools, which saves time and energy.	



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

