

ELECTRIC HOSPITAL BED - TORONCARE™ 1080



Stability meets functionality in the ToronCare™ 1080, featuring a unique four-piece stand-board foundation instead of traditional castors. This stationary electric hospital bed prioritizes safety for elderly patients and homecare settings where movement isn't required. Detachable wooden head and foot boards bring residential comfort to medical care, while powder-coated steel bed board with ventilation holes promotes air circulation. Weighing 100 kg with 2080 x 1006 mm dimensions, it serves rehabilitation centers and assisted living facilities seeking secure, non-mobile positioning.

Remote control operation manages three essential functions: back-rest elevation to 70°, leg-rest adjustment to 35°, and height range from 345-610 mm. Collapsible aluminum protective side rails incorporate user-friendly control switches for quick bedside access. The ventilated steel bed board (1960 x 850 mm) enhances patient comfort during extended rest periods, while the stand-board design eliminates rolling concerns for patients with mobility challenges or cognitive conditions requiring stationary bed placement.

Long-term care facilities gain practical advantages from this adjustable medical bed's straightforward electric controls and residential aesthetic. Standard I.V pole included, making the ToronCare™ 1080 ideal for general patient care in clinics and homes where stability outweighs mobility needs.

APPLICATIONS

- **General Patient Care** – Provides adjustable back-rest, leg-rest, and height functions for comfort and safety.
- **Hospitals and Clinics** – Suitable for wards requiring reliable, user-friendly electric beds with easy remote control.
- **Long-Term Care Facilities** – Ensures daily support with sturdy steel bed board and collapsible protective rails.
- **Rehabilitation Centers** – Assists recovery by enabling ergonomic positioning for therapy and rest.
- **Homecare Settings** – Practical for at-home patients needing professional-grade electric bed functions.
- **Elderly and Assisted Living Care** – Enhances safety with stable stand-board support instead of castors.

FEATURES

- Detachable wooden head and foot board
- Collapsible aluminum protective side rail with user-friendly control switch
- Powder-coated steel bed board with ventilated holes for increased comfort
- Bed equipped with 4 pc stand-board instead of castors, providing increased safety support for patient
- Three functions easily adjusted by remote control (back-rest, leg-rest and hi-low)

Standard Accessories (Included)

- IV Pole

TECHNICAL SPECIFICATIONS

Dimensions	
Overall dimensions (L x W x H)	2080 mm x 1006 mm x 345-610 mm (81.9 in x 39.6 in x 13.6-24.0 in)
Hi-Low	610 mm - 345 mm (24.0 in - 13.6 in)
Bed Board dimensions (L x W)	1960 mm x 850 mm (77.2 in x 33.5 in)
Bed Articulations	
Back-rest	0-70°
Leg-rest	0-35°
Weight	
Net Weight	100 Kg 220 lbs



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

