

ELECTRIC HOSPITAL BED - TORONCARE™ 1090



Simplicity drives the ToronCare™ 1090's design philosophy. This fixed-height electric hospital bed eliminates hi-low mechanisms, focusing exclusively on back-rest and leg-rest articulation through straightforward remote control operation. At just 80 kg, it's the lightest in the ToronCare electric range, making installation effortless for homecare environments and smaller clinical facilities. The fixed 380 mm height provides consistent accessibility, while detachable wooden boards deliver residential warmth to medical settings.

Back-rest reaches an impressive 75° elevation with leg-rest positioning to 45°, accommodating varied therapeutic needs without height adjustment complexity. Collapsible aluminum side rails feature integrated control switches for caregiver convenience. The powder-coated steel bed board spans 1960 x 950 mm with strategic ventilation holes that promote airflow during extended use. Overall dimensions of 2100 x 1135 mm suit standard room layouts in assisted living communities and outpatient recovery areas.

Long-term care facilities appreciate this adjustable medical bed's low-maintenance profile and cost-effective operation. Standard I.V pole included. The ToronCare™ 1090 proves that professional-grade patient support doesn't require complex mechanisms, delivering essential electric functionality for elderly care and post-surgical rehabilitation with dependable simplicity.

APPLICATIONS

- **General Hospital Wards** – Provides essential electric adjustments for patient comfort and routine care.
- **Post-Surgical Recovery** – Allows precise back-rest and leg-rest positioning to aid healing and rehabilitation.
- **Long-Term Care Facilities** – Offers durable design with collapsible side rails and ventilated steel bed board for daily use.
- **Homecare Environments** – Suitable for patients requiring professional-grade support with simple electric controls.
- **Assisted Living and Elderly Care** – Ensures comfort and safety in extended care settings with reliable, low-maintenance features.
- **Clinical and Outpatient Facilities** – Practical for patient rest and recovery in smaller healthcare centers.

FEATURES

- Detachable wooden head and foot board
- Collapsible aluminum protective side rail with user-friendly control switch
- Powder-coated steel bed board with ventilated holes for increased comfort
- Back-rest and leg-rest are easily adjusted by remote control

Standard Accessories (Included)

- IV Pole

TECHNICAL SPECIFICATIONS

Dimensions	
Overall dimensions (L x W x H)	2100 mm x 1135 mm x 380 mm (82.7 in x 44.7 in x 15.0 in)
Fixed height	380 mm (15.0 in)
Bed Board dimensions (L x W)	1960 mm x 950 mm (77.1 in x 37.4 in)
Bed Articulations	
Back-rest	0-75°
Leg-rest	0-45°
Weight	

Net Weight	80 Kg (176 lbs)
------------	-----------------



TORONTECH

Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com



USA | Canada | EU | UAE | GCC | KSA | India | Africa | Latin America | Asia

