



EXHALATION VALVE AIR LEAKPROOFNESS TESTER

Exhalation Valve Air Leakproofness Tester ToronFM™ EVALT is a high-precision device designed for accurate air leak proofness testing. It features a differential pressure gauge, a user-friendly pressure adjustment system, and eliminates timing inaccuracies common with manual testing methods.

APPLICATIONS

Standards

- 1. EN 149
- 2. NIOSH - 42 CFR Part 84


FEATURES

- 1. Advanced air leakproofness testing device designed to eliminate timing inaccuracies caused by manual stopwatch methods.
- 2. High-precision, differential pressure gauge for accurate results.
- 3. User-friendly and efficient pressure adjustment system for seamless operation.

TECHNICAL SPECIFICATIONS

Model	ToronFM™ EVALT
Bumper Capacity	≥5L

Model	ToronFM™ EVALT
Micro Pressure Difference Gauge Range	0~2000Pa, Accuracy 1Pa
Flowmeter Range	0~1000ml/min, Accuracy 3%
Pumping Rate	2L/min
Voltage	AC220V 50Hz 150W (110V is also available)
Dimensions	500x250x450mm



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

