



TORONTECH



PRESET TYPE BLOOD SPECIMEN COLLECTOR (BLOOD COLLECTING PERFORMANCE TESTER)

Preset Type Blood Specimen Collector (Blood Collecting Performance Tester) ToronFM™ BCPT is designed for single-use human arterial blood collectors. It ensures proper piston functionality before blood collection, checks for synthetic blood leakage at the end of the collection process, and verifies the system's airtightness after the procedure.


APPLICATIONS

Standard

- 1. YY 0612 (Section 4.4) – Specifies the particle filtration efficiency test requirements for medical protective masks, defining how effectively the mask filters airborne particles under standardized test conditions.

TECHNICAL SPECIFICATIONS

Model	ToronFM™ BCPT
Pressure of Hydrostatic Head	700mm (liquid column)
Inner Diameter of Channel	≥2.5mm (liquid reservoir connects with channel)
Water Cycle System	(37±1) °C
Pressure of Liquid Reservoir	50kPa
Timer	0-99.99s/m/h, Accuracy: ±0.01
Standards	YY 0612 (4.4)



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

