



WASTE GAS RECYCLER

Waste Gas Recycler efficiently captures discarded doses during delivered dose uniformity testing for inhaled MDI, spray, and DPI products. This compact desktop device ensures safe and convenient collection, safeguarding laboratory staff from exposure to waste drugs.


APPLICATIONS

Standards

- United States Pharmacopeia (USP) <601> - Aerosols, Nasal Sprays, Metered-Dose Inhalers, and Dry Powder Inhalers
- European Pharmacopoeia (Ph. Eur.) 2.9.18 - Preparations for Inhalation: Aerodynamic Assessment of Fine Particles
- Chinese Pharmacopoeia (ChP) 2020 - Inhalation Product Testing
- Laboratory Safety Standards for Occupational Exposure Prevention

FEATURES

1. Filter: capsule type.
2. Filter material: fiberglass.
3. Filtration efficiency: 99% for particles above 0.2 microns.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

