



TORONTECH



SURGICAL TABLE - FULLY ELECTRIC TORONSURGE™ -1000A SERIES

The ToronSurge™-1000A is a state-of-the-art fully electric surgical table, available in three variations to meet diverse surgical needs:

1. S – Standard
2. E – Economical
3. D – Deluxe

APPLICATIONS

Powered by an electric push rod motor operating system, the ToronSurge™-1000A allows seamless adjustments to accommodate various surgical positions, including:

1. Table height adjustment
2. Forward and backward tilt
3. Left and right tilt
4. Back plate folding
5. Tabletop translation

This multifunctional electric operating table is ideal for a wide range of surgical procedures, including abdominal surgery, obstetrics, gynecology, ENT, urology, anorectal, orthopedics, and more.

[Surgical Table - Fully Electric ToronSurge™ -1000A Series](#)

[Fully Electric Surgical Table - ToronSurge™ -1000A Series](#)

Fully Electric Surgical Table - ToronSurge™ -1000A Series

The Surgical Table - Fully Electric ToronSurge™ -1000A Series is a versatile, high-performance solution designed to meet the evolving needs of modern surgical environments. With X-ray compatibility, dual control options, a powerful battery, and an integrated kidney bridge, this table ensures efficiency, safety, and reliability across a wide range of medical procedures.

FEATURES

Key Features of the ToronSurge™ -1000A

1. X-ray Compatibility for Enhanced Imaging

Designed with a PFCC tabletop, the ToronSurge™ -1000A Series allows for X-ray scanning during procedures. The ToronSurge™ -1000A series offers a tabletop translation of over 300mm, providing excellent C-arm imaging access and compatibility with X-ray film boxes.

2. Dual Control System for Reliability

Equipped with a hand controller and an optional control panel, this system ensures enhanced surgical safety and flexibility, providing redundant control options for uninterrupted operation.

3. High-Performance Rechargeable Battery

The built-in rechargeable battery supports up to 50 operations on a single charge, ensuring reliable performance even in cases of power outages. Additionally, AC power supply compatibility ensures continuous electrical energy for maximum safety and efficiency.

4. Built-in Kidney Bridge for Specialized Procedures

Featuring an integrated lumbar bridge, the table provides optimal positioning for bile duct and kidney surgeries, improving surgical precision and ease of access.

TECHNICAL SPECIFICATIONS

Specification	Standard Version	Economic Version	Deluxe Version
Length & Width	2070mm × 550mm		
Elevation (Up/Down)	1000mm / 700mm		
Head Plate (Up/Down)	45° / 90°		
Back Plate (Up/Down)	75° / 20°		
Leg Plate (Up/Down/Outward)	15° / 90° / 90°		
Trendelenburg/Reverse Trendelenburg	25° / 25°		
Lateral Tilt (Left/Right)	15° / 15°		
Kidney Bridge Elevation	≥110mm		
Horizontal Sliding	300mm		
Flex/Reflex	Combination Operation		
X-ray Compatibility	Optional		
Control Panel	Standard		
Electro-motor System	Electric Motor		
Voltage	220V / 110V		

Specification	Standard Version	Economic Version	Deluxe Version
Frequency	50Hz / 60Hz		
Power Capacity	1.0 KW		
Battery	Yes		
Mattress	Memory Mattress		
Main Material	304 Stainless Steel		
Maximum Load Capacity	200 KG		
Warranty	1 Year		
Double Control System	Optional	Optional	Yes
Leg Plate Adjustment	Manual	Manual	Pneumatic

Standard Accessories

Accessories for Enhanced Functionality



Anaesthesia Screen

Provides a privacy barrier during procedures.



Fixing Clamp

Securely attaches accessories to the table rail.



Stomach Support Rod

Supports abdominal positioning during surgery.




IV Stand

Holds IV bags for fluid administration.



Fixing Clamp

Securely attaches accessories to the table rail.



Knee Support Rod

Provides support and stability for the knee.




Arm Rest (PU)

Offers comfortable support for the arm.




Fixing Clamp

Securely attaches accessories to the table rail.



Seamless Memory Pad

Conforms to the body for enhanced comfort.




Arm Rest (Padded)

Offers comfortable support for the arm.



Fixing Clamp

Securely attaches accessories to the table rail.



Memory Pad

Conforms to the body for enhanced comfort.




Arm Rest (Padded)

Offers comfortable support for the arm.




Fixing Clamp

Securely attaches accessories to the table rail.



PU Pad

Durable surface for added support.



Head Ring

Provides head support and stability.



Body Belt (Buckle Type)

Secures the body for safety and stability.



Waste Basin

For convenient collection of waste.



Lateral Support Pad

Provides lateral support and positioning stability.




Body Belt

Secures the body for safety and stability.



Operation Support

Enhances support and positioning during procedures.



Prone Position Pad

Supports the body in prone position for procedures.



Leg Belt

Secures the legs to help maintain proper positioning.




Touchscreen Hand Remote

Intuitive touchscreen control for table adjustments.



Stomach Support Pad

Supports the abdomen to improve patient comfort.



Arm Belt

Secures the arms to reduce movement and enhance safety.



Button Hand Remote


Easy-to-use button control for table positioning.

Accessories for Surgical Table

**Leg Support
(Padded)**



Wrist Belt



**Hand
Remote**



**Leg Support
(PU)**



Foot Plate



**Lateral
Arm Support**



**External
Kidney Bridge**



**Shoulder
Support**



**Instrument
Tray**



Arm Surgery
Stand



Body
Support



Orthopedic
Traction 10 05



Clamp



Body
Support




Orthopedic
Traction 10 06



Clamp




Head Rest




Three Point
Skull Clamp





Accessory Name	Quantity
Anesthesia Screen	1 piece
Body Support	1 pair
Arm Support	1 pair
Shoulder Support	1 pair
Leg Support	1 pair

Accessory Name	Quantity
Foot Support	1 pair
Kidney Bridge Handle	1 piece
Mattress	1 set
Fixing Clamps	8 pieces
Long Fixing Clamps	1 pair
Hand Remote	1 piece
Power Line	1 piece



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

