



VISCOMETERS

Viscometers for Accurate Fluid Viscosity Testing

Achieve absolute confidence in your fluid characterization. The ToronVisc™ series of advanced viscosity meters provides the verifiable data you need to ensure product consistency, optimize formulations, and protect your return on investment by eliminating out of spec materials. These robust instruments are engineered for demanding lab and production environments where precision is not optional.

Our ToronVisc™ instruments deliver the accuracy and repeatability you expect from premium systems. They are fully compatible with established testing methods, including those for Brookfield-style viscometers, providing a direct and cost-effective pathway to achieving industry standard results. When your process demands reliable rheological data for Newtonian or non-Newtonian fluids, our viscosity measurement tools provide an ideal balance of performance and value.

Selecting the Right Viscosity Meter for Your Application

Torontech provides a complete portfolio of viscosity measurement tools designed for specific testing requirements.

- **Rotational Viscometers:** The versatile workhorse for multipurpose QC and general materials testing.
- **Cone and Plate Viscometers:** The precision instrument for small sample volumes and applications requiring absolute shear rate data.
- **Digital and Programmable Models:** Automate complex test procedures and enhance data management with intuitive touchscreen controls.

Engineered for Your Industry's Demands

Our line of ToronVisc™ viscosity meters is trusted by industry leaders for critical measurements in:

- Paints, Coatings, and Inks

- Cosmetics and Personal Care
- Pharmaceuticals and Biotech
- Food and Beverage
- Petrochemical Products and Lubricants
- Adhesives and Construction Chemicals
- Polymer and Chemical Formulations

Strategic Advantages and Key Features

Every ToronVisc™ instrument is designed with features that deliver tangible benefits to your operation.

- **Expansive Measurement Capability:** A wide viscosity range allows you to confidently measure everything from thin solvents to thick pastes.
- **Intuitive Digital Interface:** Reduce operator learning curves and minimize the potential for manual errors.
- **Versatile Spindle Sets:** Adapt to diverse fluid types and test requirements with included interchangeable spindles for maximum flexibility.
- **Seamless Data Integration:** Export results via USB or RS232 directly to your LIMS or statistical process control software for enhanced traceability.
- **Precise Temperature Control:** Add optional temperature units to eliminate a critical variable and guarantee repeatable, comparable results.
- **Guaranteed Standards Compliance:** Fully conform to ASTM D2196, ISO 2555, ISO 1652 for validated, auditable reporting.

As a global supplier leveraging our innovative technologies, Torontech engineered the ToronVisc™ line to challenge the industry's price to performance benchmarks. We believe precise, reliable viscosity measurement should be accessible, not exclusive. When your process control and quality assurance demand consistent data, trust ToronVisc™ viscosity meters to protect your product and your budget.

Footnote/Disclaimer:

Brookfield is a registered trademark of AMETEK, Inc. All other product and company names are trademarks™ or registered®

trademarks of their respective holders. Use of them does not imply any affiliation with or endorsement by them.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

