


INTELLIGENT ALL-IN-ONE VISCOMETER - TORONVISC™ T1



Advanced Touchscreen Viscometer with Integrated Temperature Control

ToronVisc™ T1 Intelligent All-in-One Viscometer is an innovative solution for high-precision viscosity measurement across a wide range of applications. This fully integrated system features a 7-inch touchscreen interface, built-in temperature control, and automatic electric lift for seamless operation and real-time rheological monitoring.

Designed to simplify complex workflows in R&D and quality control labs, the ToronVisc™ T1 combines intelligent software, precise thermal regulation, and high measurement stability into one compact, easy-to-use unit. Whether you're measuring low, medium, or high viscosity fluids, this all-in-one viscometer ensures reliable performance and consistent data acquisition.

APPLICATIONS

Applications

- Ideal for low, medium, and high viscosity fluid testing
- Suited for industries such as cosmetics, food & beverage, pharmaceuticals, petroleum, and chemicals
- Enables temperature-dependent viscosity profiling with integrated heating/cooling

Optional Accessories

- SNB Viscosity Analysis Software
- Certified Viscosity Standard Oils

Disclaimer: Brookfield® is a registered trademark of AMETEK, Inc. Torontech's ToronVisc™ series is designed to be fully compatible in performance and functionality, but is not affiliated with or endorsed by AMETEK® or Brookfield®.

FEATURES

- 7" Touchscreen Display - Intuitive, rich interface for easy data input and monitoring
- Electric Lift & Auto Positioning - Ensures spindle consistency and operational safety
- Shear Rate and Stress Displayed - Ideal for comprehensive rheological analysis
- Integrated Temperature Control - Single-point temperature setting from -1°C to 150°C
- Data Management - External storage support (single point, continuous, timed)
- Calibration Functionality - Password-protected temperature and viscosity calibration
- Connectivity - USB port and Gigabit Ethernet for high-speed data transfer
- Optional SNB Software - Record and export viscosity data for in-depth analysis
- Stable Power Supply - Operates on 100V-240V with strong anti-interference performance

TECHNICAL SPECIFICATIONS

Parameter	Specification
Model	ToronVisc™ T1
Viscosity Range	L: 5-0.33M mPa·s R: 50-3.3M mPa·s HA: 100-6.6M mPa·s HB: 400-26.4M mPa·s
Temperature Range	-1°C to 150°C
Rotational Speed	0.3-100 RPM
Spindles	SC4-21, 27, 28, 29 (standard) Optional: SC4-14, 15, 16, 18, 25, 31, 34
Accuracy	±1.0% (Full Scale)

Parameter	Specification
Repeatability	±0.2% (Full Scale)
Voltage	AC100-240V, 50Hz/60Hz Output: DC12V/1.5A



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

