



TORONTECH



OPTIONAL ACCESSORIES FOR TORONVISC™ VISCOMETERS

Enhance the performance and versatility of your ToronVisc™ viscometers with a full range of high-quality optional accessories. These add-ons are designed to support advanced testing scenarios, ensure precise data collection, and improve user efficiency across various viscosity measurement applications.

APPLICATIONS

Optional Accessories for ToronVisc™ Viscometers

Ultra Low Viscosity Adapter (ULA)

Optional Accessories for ToronVisc™ Viscometers

Small Sample Adapter (SSA)

Optional Accessories for ToronVisc™ Viscometers

Micro Printer

Optional Accessories for ToronVisc™ Viscometers

RTD and PT100 Temperature Probes

Ultra Low Viscosity Adapter (ULA)

- Ideal for measuring ultra-low viscosity fluids (as low as 1 mPa.s)
- Sample volume: 30 mL
- Simple attachment to a standard viscometer
- Compatible with advanced models (e.g., NTV-P3) for shear rate testing

Small Sample Adapter (SSA)

- Designed for sample volumes between 5–18 mL
- Easy to install and clean
- Delivers accurate shear rate and stress data for detailed product analysis
- Made of 304 or 316 stainless steel for durability and chemical resistance

Thermostatic Water Bath

- Temperature range: -5°C to +95°C
- Cooling system: rotary compressor
- Temperature stability: $\pm 0.05^{\circ}\text{C}$; Accuracy: 0.1°C
- Accommodates one or more sample containers

SNB Software

- Automatic data collection and viscosity curve creation
- Compare historical data (up to 10 sets)
- Export records to Excel format
- Ideal for quality control and R&D labs

DV Software

- Enables automatic data collection and measurement
- Generates viscosity curves and printable graphs
- Compare up to 10 historical datasets
- Offline work support and Excel export capability

Micro Printer

- Simple, user-friendly operation
- Factory-ready setup
- Prints real-time values for non-Newtonian fluid measurements

RTD and PT100 Temperature Probes

- RTD temperature range: 0°C to 100°C
- PT100 temperature range: -50°C to 200°C
- Ensures precise real-time temperature readings

Standard Oils

- Used to verify accuracy and calibrate viscometers
- Available in a broad range of viscosities:

- 5 Series, 200 Series, 500 Series, 1000 Series
- 3000 Series, 6000 Series, 10000 Series
- 30000 Series, 60000 Series, 100000 Series

Whether you're performing routine quality checks or in-depth rheological research, these accessories provide the flexibility and precision needed for optimized viscosity measurement with the ToronVisc™ series.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

