



WATER PURIFICATION EQUIPMENT

Water Purification Equipment / Water Purification System from Torontech offers tap water inlet, integrating various advanced water purification technology and process control methods, our water purification systems can completely meet, even exceed, the pure water standard of ASTM, NCCLS Class I etc. The ultra pure water resistivity is assured to reach 18.2 M Ω -cm.

It can meet the need of advanced lab water use for a different purpose. It can be widely used in the field of chemistry, biology, pharmacy, medicine, microelectronics, semiconductor etc, meeting the application requirements of spectral analysis, chromatographic analysis, cell culture, protein purification and molecular biology.

Torontech offers Water Purification Equipment / Water Purifiers models: TT Alpha,TT Beta,TT Delta,TT Gamma,TT Kappa,TT Omega,TT Sigma &TT Zeta

Advantages

- The super-large LCD display, display the system running state and various parameters intuitively
- Automatic microcomputer controlling system, multi-menu operating, animation mode display
- Fault automatically detect, automatic diagnosis
- Water quality over standard alarm, no water alarm, consumables end alarm function
- Consumables residual life shows, inform the user to replace consumables timely
- 3 road on-line water quality monitoring, monitor the inlet, RO and Ultra pure water's quality timely
- Built-in system disinfection procedure, realize the disinfection of ultra pure water's tube
- Embedded sterilization program to achieve full pipeline disinfection

Torontech offers Water Purification Equipment / Water Purifiers models: TT Alpha,TT Beta,TT Delta,TT Gamma,TT Kappa,TT Omega,TT Sigma &TT Zeta



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

