



BIRTHING BEDS

Exploring Birthing Bed, Gynecological Examination Tables and Obstetric Beds ToronComfort™ Series from Torontech

Manual, hydraulic, and electric birthing beds for labor, delivery, and recovery, built to IEC 60601-2-52.

A birthing bed supports a mother through labor, delivery, gynecological examination, and recovery on a single platform, helping obstetrics teams manage each stage without moving the patient. Torontech offers manual, hydraulic, and electric models built to IEC 60601-2-52, the safety standard for medical beds. The range spans simple examination beds to full labor-delivery-recovery (LDR) systems.

Maternity hospitals, birthing centers, and gynecology clinics choose a configuration based on acuity, room size, and workflow. A manual bed suits routine exams and lower-acuity births, while an electric LDR bed handles a busy delivery suite. This page compares the options so you can match the bed to your department.

Why Choose Torontech for Your Birthing Bed

Torontech supplies testing and measurement equipment to laboratories and manufacturers across North America and worldwide, backed by responsive technical support and calibration service. The birthing bed range is built for reliable, standards-driven results, so you can standardize on one supplier as your testing needs grow. Our team helps you match the right configuration to your samples, applicable standards, and throughput, then supports you through installation and beyond.

Request a Quote

Tell us your delivery volume, the positions you need, and whether you want a manual, hydraulic, or electric bed, and we will recommend the right birthing bed. [Contact the Torontech team](#) for pricing and lead times.

APPLICATIONS

Applications

Birthing beds serve the full obstetric and gynecological care pathway across hospitals and clinics.

- Hospital Labor and Delivery: carries a mother through labor, birth, and early recovery, reducing transfers during active labor
- Gynecology Clinics and Exam Rooms: provides a stable, adjustable surface for pelvic exams, prenatal checks, and minor procedures
- Maternity Wards and Postpartum Care: supports comfortable positioning for rest, monitoring, and breastfeeding after birth
- Birthing Centers: offers a clinical feature set in a patient-friendly form for lower-acuity, natural births
- Obstetric Emergencies: rapid flat or Trendelenburg positioning supports urgent care and caesarean preparation

TECHNICAL SPECIFICATIONS

Birthing Bed Range

Torontech birthing beds fall into two families: the manual and hydraulic ToronComfort™ 1000 Series and the electric ToronComfort™ 1010. The table below compares the four models so you can match drive type and positioning to your clinical needs.

Model	Type	Adjustment	Best For
-------	------	------------	----------

ToronComfort™ 1000-A	Hydraulic	Height 715-965 mm, backrest to 55°, Trendelenburg tilt	Delivery suites needing height and clinical tilt without electric power
ToronComfort™ 1000-B	Manual, fixed height (narrow)	Adjustable backrest, removable leg plate	Gynecology examination where a narrower 620 mm frame fits the room better
ToronComfort™ 1000-C	Manual, fixed height	Adjustable backrest, removable leg plate	General labor and exam rooms needing a dependable exam and delivery bed
ToronComfort™ 1010	Electric	Height 485-860 mm, 15+ birth positions, mechanical CPR	Maternity hospitals running a full LDR workflow in one room

How to Choose the Right Birthing Bed

- For routine gynecology exams and lower-acuity births: choose a manual model like the ToronComfort™ 1000-C or ToronComfort™ 1000-B for a dependable, lower-cost surface.
- For tight exam rooms: choose the narrower ToronComfort™ 1000-B to fit a 620 mm frame where space is limited.
- For delivery suites that need height and tilt without power: choose the hydraulic ToronComfort™ 1000-A for a gas-spring backrest, Trendelenburg positioning, and mobility.
- For a full labor-delivery-recovery workflow: choose the electric ToronComfort™ 1010 for powered positioning, more than fifteen birth positions, and mechanical CPR.
- For mobility between rooms: choose the ToronComfort™ 1000-A or ToronComfort™ 1010, both fitted with central-brake casters.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

