



SURGICAL LIGHT - SURGICAL LAMP

ToronLight™ Series & Surgical Microscopes — Illuminate Every Surgical Moment

Surgical Lights — also known as surgical lamps or operating-theater lights — and Surgical Microscopes together form the visualization foundation of every modern operating room. Without bright, shadow-free, cool-spectrum illumination and the precision optics that magnify the smallest anatomical structures, even the most skilled surgical team cannot consistently deliver world-class results. From routine procedures in general surgery to the most delicate microsurgical interventions in ophthalmology, ENT, and neurosurgery, surgical lighting and visualization equipment directly determines what the surgeon can see — and therefore what the surgeon can do.

Torontech proudly presents the ToronLight™ 1000 Series surgical lights and a comprehensive range of Surgical Microscopes engineered to meet the diverse visualization needs of modern surgical practice. The ToronLight™ 1000 is a state-of-the-art LED surgical light featuring zero-heat cool-light technology, shadow-free reflector design, large depth of field, and a slim compact light head mounted on a 7-axis mechanical arm for exceptional positioning flexibility. Complementing it, Torontech's Surgical Microscopes deliver high-magnification apochromatic optics, multi-step or motorized zoom, ergonomic binocular surgeon heads with optional assistant heads, integrated visualization (CCD cameras, video recording, beam splitters), and specialized configurations for ophthalmic, ENT, neurosurgical, dental, and plastic-reconstructive applications. Together, the two product families equip every operating room — from the largest tertiary hospital to the most specialized surgical center — with the illumination and visualization performance that modern surgery demands.

See Every Detail. Operate with Confidence.

In every operating room, what the surgeon can see defines what the surgeon can do. The right surgical light illuminates the deepest cavity without shadows, without heat, without strain — for hours on end. The right surgical microscope magnifies the smallest blood vessel, the finest nerve, the tiniest structure with the clarity and stability the procedure demands. The hospitals and surgical centers that consistently deliver world-class outcomes are the ones equipped with surgical

visualization built for precision, comfort, and reliability — and the Torontech ToronLight™ surgical lights and Surgical Microscopes give you exactly that.

Whether you are equipping a new operating room, upgrading an outdated lighting system, or specifying a high-end microscope for an ophthalmology or neurosurgery suite, Torontech has the visualization platform engineered to meet your exact needs. Equip your operating room with the surgical lighting and microscope solutions trusted by surgeons, OR managers, and biomedical engineers worldwide.

APPLICATIONS

Applications in Surgical Specialties

- **General & Abdominal Surgery** — Bright, shadow-free OR illumination for laparotomy, hernia, GI, and broad-scope general procedures.
- **Ophthalmology / Eye Surgery** — High-magnification surgical microscopes with retina-protection illumination, red-reflex enhancement, and X-Y motorized centering for cataract, vitreoretinal, and refractive surgery.
- **ENT (Otolaryngology)** — Compact apochromatic surgical microscopes for ear, nose, and throat surgery, otologic procedures, and microlaryngoscopy.
- **Neurosurgery** — Stable, vibration-free, high-magnification visualization combined with deep, shadow-free LED illumination for cranial, spinal, and neurovascular interventions.
- **Plastic, Reconstructive & Hand Surgery** — Microscope-assisted microvascular and microneural surgery, free-flap reconstruction, and fine anastomosis under shadow-free cool-light OR illumination.
- **Dental & Oral Surgery** — Operating dental microscopes for endodontics, oral surgery, and implantology requiring high-magnification visualization of root canals and small structures.
- **Orthopedic, Cardiovascular & OBGYN** — ToronLight™ 1000 illumination for all major-specialty operating rooms, day-surgery centers, emergency theaters, and ambulatory surgical suites.

FEATURES

Why Choose Torontech for Your Surgical-Visualization Needs

With more than two decades of experience supplying advanced medical equipment to hospitals and surgical centers across more than 60 countries, Torontech delivers premium surgical-illumination and visualization performance backed by global expertise.

Engineered for the Surgical Environment

The ToronLight™ 1000 Series is designed and manufactured to the highest industry standards. Zero-heat LED emission keeps OR staff comfortable during long procedures, shadow-free reflector design ensures consistent visualization of deep cavities, and the 7-axis mechanical arm gives surgeons full positioning freedom from any angle around the operating table.

Integrated Lighting + Visualization

By offering both surgical lights and surgical microscopes under one trusted brand, Torontech equips surgical teams with a fully coordinated visualization ecosystem — from broad operating-field illumination down to the highest-magnification microsurgical detail — simplifying procurement, service, and training.

Energy Efficiency & Long Life

ToronLight™ LEDs consume considerably less power than halogen alternatives while delivering long-life service that virtually eliminates the disruption of bulb changes mid-shift. Combined with cool-light operation, the result is lower energy costs, reduced HVAC load on the OR, and improved comfort for the entire surgical team.

Cost-Effective Value & Global Support

Torontech delivers premium surgical-light and microscope performance at a transparent, competitive price point — without the inflated cost structures of legacy brands. With offices and authorized partners across the USA, Canada, UAE, KSA, GCC, EU, India, APAC, Africa, and Latin America, expert installation, surgical-team training, and ongoing service are always within reach.

TECHNICAL SPECIFICATIONS

Product Specification Table

Specification	ToronLight™ 1000 Surgical Light	Surgical Microscopes
Product Type	LED Surgical Light (operating-theater illumination)	Operating microscope for high-magnification microsurgery
Primary Function	Shadow-free, cool-light illumination of the surgical field	High-magnification visualization of fine anatomical structures
Light Source / Optics	Long-life LED — generates no heat; energy-efficient cool light	Apochromatic optical system with manual or motorized zoom; high-grade imported glass for minimal aberration
Magnification	Not applicable	Multi-step magnification (typical surgical microscope range 4×–40×)

Shadow / Image Control	Reflectors effectively remove shadows; large depth of field	Beam splitter / illumination control; retina-protection and red-reflex enhancement on ophthalmic models
Positioning / Movement	Mechanical arm with 7 axes of rotation — exceptionally large movement envelope; light head positionable from any angle	Multi-axis articulating arm; foot-controlled or motorized X-Y centering and focus on advanced models
Heat Output	Zero heat — LEDs emit cool light, keeping OR staff and patient comfortable during long surgeries	Low — LED illumination on modern models
Energy Efficiency	Significantly reduced energy consumption vs. halogen / incandescent	Energy-efficient LED illumination
Visualization Add-Ons	Single, dual, or triple-dome configurations; HD video camera options	Binocular surgeon head + binocular or monocular assistant head; CCD camera adapter + video recording options
Sterility / Hygiene	Sterilizable removable handles; smooth, easy-clean light head surfaces	Sterilizable handle covers; smooth surfaces for OR disinfection protocols
Mounting Options	Ceiling-mounted (standard); wall and mobile-stand variants available	Floor, ceiling, or wall-mounted configurations
Surgical Specialties	All OR & procedure rooms — general, OBGYN, ortho, ENT, cardio, neuro, day-surgery	Ophthalmology, ENT, neurosurgery, dental, plastic & reconstructive, hand surgery
Compliance	Designed & manufactured to highest industry standards (IEC-aligned)	International medical-device manufacturing & optical standards



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

