



TORONTECH

MEDICAL GLOVES WATERTIGHTNESS TESTER



Medical Gloves Watertightness Tester | ISO 10282, EN 455-1, ASTM D5151

Medical Gloves Watertightness Tester – Model ToronMed™ MGWT01 is designed to evaluate the watertightness of medical examination gloves and protective gloves used in healthcare, laboratory, chemical, and microbiological environments.

By accurately simulating real-world use conditions, this tester ensures that gloves provide reliable liquid barrier protection against contamination, infection, and hazardous exposure. It is an essential tool for glove manufacturers, regulatory laboratories, and healthcare product quality control.

Standards

The Medical Gloves Watertightness Tester – Model ToronMed™ MGWT01 is designed in compliance with relevant international standards for medical and protective gloves, including:

- ISO 10282 – Single-use sterile rubber surgical gloves – determination of watertightness.
- EN 455-1 – Medical gloves for single use – requirements and testing for freedom from holes.
- ASTM D5151 – Standard test method for detection of holes in medical gloves.

By meeting these standards, the ToronMed™ MGWT01 ensures glove manufacturers and testing laboratories can demonstrate global compliance, safety, and reliability of protective gloves.

The Medical Gloves Watertightness Tester – Model ToronMed™ MGWT01 provides a cost-effective, reliable, and standards-compliant solution for testing the liquid barrier integrity of gloves, ensuring they deliver maximum protection and quality assurance for both medical professionals and patients.

APPLICATIONS

The Medical Gloves Watertightness Tester – Model ToronMed™ MGWT01 is suitable for:

- Medical examination gloves used in hospitals, clinics, and laboratories.
- Protective gloves for handling contaminated materials.
- Chemical-resistant gloves used in pharmaceutical and laboratory settings.
- Microbial barrier gloves designed for infection prevention and biosafety.

FEATURES

- Specialized glove test fixtures for high-efficiency testing of multiple glove types.
- Transparent hollow mandrels (≥ 60 mm diameter) designed to accommodate up to 1000 mL of water for watertightness evaluation.
- Sturdy holding device to ensure stability and accuracy during testing.
- Graduated cylinder (≥ 1000 mL capacity) for precise water measurement.
- Integrated timer to ensure accurate test cycles.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

