

CHEMICAL PROPERTIES & TITRATION ANALYZERS



Chemical Properties & Titration Analyzers

Torontech's Chemical Properties & Titration Analyzers category brings together a comprehensive lineup of precision analytical instruments engineered to quantify the chemical composition, stability, and contamination levels of petroleum products, lubricants, fuels, and industrial fluids. These instruments form the analytical backbone of modern petrochemical refineries, lube oil blending facilities, fuel terminals, and quality control laboratories worldwide.

Chemical characterization is the foundation of product quality assurance in the petroleum and petrochemical industry. Whether verifying the acid neutralization capacity of a fresh engine oil, detecting water ingress in a transformer fluid, determining the salt content of crude oil before refinery processing, or confirming the sulfur and chlorine levels of a specialty chemical, the integrity of the titration and chemical analysis result directly determines product acceptance, safety, and regulatory compliance.

Torontech's portfolio in this category spans potentiometric and coulometric titration systems, Karl Fischer water analyzers, elemental analyzers for sulfur and chlorine, acid and base number testers, oxidation stability instruments, and multifunctional chemical property analyzers — all designed and validated against internationally recognized ASTM, ISO, IP, and IEC test methods.

The Role of Titration Analysis in Petrochemical Quality Control

In the petroleum industry, chemical properties testing is not optional — it is mandatory at nearly every stage of the production and distribution chain. Crude oil must be characterized for salt, water, and sulfur content before entering the refinery desalter and distillation columns. Refined products must meet strict acid number, sulfur, and contaminant specifications before release to market. Lubricating oils must be tested for base number depletion during in-service monitoring to determine change intervals. Transformer oils require water content and acidity analysis to predict dielectric breakdown risk.

Torontech’s Chemical Properties & Titration Analyzers are purpose-built to deliver the sensitivity, selectivity, and reproducibility that these mission-critical quality decisions demand — whether operating in a centralized refinery laboratory, a portable field-testing unit, or a lube oil formulation R&D center.

APPLICATIONS

Applications — Petrochemical & Industrial Focus

Torontech’s Chemical Properties & Titration Analyzers are deployed across the full breadth of the petroleum and petrochemical value chain — from upstream exploration and crude processing through refining, blending, finished product quality control, and in-service fluid monitoring. The following application areas represent the core deployment environments for this category.

Industry Sector	Key Chemical Analysis Applications
Crude Oil & Upstream	Salt content of crude oil (ASTM D3230) for desalter control; water & sediment analysis; acidity screening for pipeline corrosion risk assessment; mercaptan sulfur for sour crude characterization.
Petroleum Refining	Acid number monitoring of process streams; total sulfur in naphtha, diesel, jet fuel, and heavy fuel oils (ASTM D5453, D7039); chlorine in feedstocks to protect catalyst beds; water content in refined products.
Lubricant Manufacturing & Blending	Total base number (TBN) of engine oils, marine diesels, and metalworking fluids (ASTM D2896, D4739); acid number of fresh and in-service oils (ASTM D664); oxidation stability of lubricant base stocks (ASTM D2272); hydroxyl and amine values of additive packages.
Fuel Quality & Distribution	Acid number of aviation turbine fuels and biodiesel (ASTM D664, D974); water content of jet fuel, gasoline, and diesel (Karl Fischer, ASTM D6304); mercaptan sulfur of natural gas liquids (ASTM D3227); chlorine in finished fuel blends.

Industry Sector	Key Chemical Analysis Applications
Transformer & Electrical Oils	Neutralization number / acidity of transformer, turbine, and hydraulic oils (ASTM D974, IEC 60296); Karl Fischer water content for dielectric strength prediction; oxidation stability for service life estimation.
Petrochemical Intermediates	Acid value, hydroxyl value, and amine value of polyols, fatty acids, and chemical intermediates; chlorine and sulfur content of raw materials and specialty chemicals; saponification number of ester-based products.
In-Service Oil Monitoring	Trending of acid number increase and TBN depletion in engine, turbine, and compressor oils for condition-based maintenance programs; water ingress detection; oxidation degradation monitoring.
Marine & Bunker Fuel	Acidity and TBN of marine cylinder lubricants; sulfur content of bunker fuel oils for IMO 2020 compliance (ASTM D5453); salt and water contamination screening in fuel bunkering operations.

Petrochemical Laboratory Workflows Supported

Torontech's instruments are designed to integrate seamlessly into the following laboratory workflow environments:

QC Lab Incoming & Outgoing Product QC	High-throughput acid number, water content, and TBN testing for finished lubricants, fuels, and petrochemical products at production release checkpoints. Designed for rapid cycle times and minimal sample preparation.
Refinery Refinery Process Lab	Real-time monitoring of chemical properties along refinery process streams — salt in desalted crude, chlorine in naphtha, acidity in distillate fractions — to enable timely process corrections and prevent downstream equipment damage.
R&D Lubricant R&D & Formulation	Precise potentiometric and coulometric titration for additive package characterization, base oil screening, and blended formulation validation against ASTM and OEM specifications.

Field On-Site & Field Testing	Compact and portable chemical analyzers for in-service oil sampling and rapid titration in power generation facilities, marine vessels, and industrial plant maintenance programs.
Env Environmental & Regulatory Testing	Sulfur, chlorine, and nitrogen analysis of petroleum products and feedstocks for environmental emissions compliance and regulatory reporting under EPA, REACH, and international standards.

FEATURES

Why Choose Torontech for Chemical Properties & Titration Analysis

Torontech is a specialized analytical instrument provider with deep domain expertise in petroleum and petrochemical testing. When you choose Torontech for your chemical properties and titration analyzer needs, you gain more than hardware — you gain a partnership with a team that understands the unique analytical challenges, regulatory landscape, and operational demands of the petroleum industry.

ASTM Full ASTM & ISO Method Compliance	Every instrument in this category is designed, calibrated, and validated against the specific ASTM, ISO, IP, and IEC test methods that govern petroleum product quality in North American and international markets — including ASTM D664, D2896, D4739, D974, D6304, D5453, D7039, D3230, D2272, and D3227. Compliance documentation and method correlation data are available for all supported standards.
HC Petroleum-Specific Application Expertise	Unlike generic laboratory instrument suppliers, Torontech's application engineers specialize in petroleum and petrochemical testing. We understand the matrix challenges of crude oil and heavy fuel oil samples, the electrode management requirements of high-TBN and high-acid matrices, and the Karl Fischer optimization needed for complex petroleum-water systems. This domain expertise is embedded in every instrument's default method library and application support documentation.

DATA Integrated Digital Data Management	All Torontech chemical analyzers generate traceable, audit-ready result records with operator ID, timestamp, sample ID, and full titration curve data. Results can be exported in standard formats compatible with LIMS, ERP, and QC data management systems — supporting 21 CFR Part 11, ISO/IEC 17025, and refinery quality management system requirements.
RANGE Complete Category Coverage in One Supplier	Torontech provides a single-source solution covering potentiometric auto-titrators, Karl Fischer volumetric and coulometric water analyzers, elemental sulfur and chlorine analyzers, oxidation stability testers, and multi-parameter chemical analyzers. This eliminates the complexity of managing multiple instrument suppliers, diverse service contracts, and incompatible calibration schedules across your chemical analysis laboratory.
SVC Application Support & Method Development	Our technical team provides instrument commissioning, operator training, IQ/OQ/PQ qualification support, and ongoing application assistance for method development and troubleshooting. For petroleum laboratories implementing new ASTM methods or transitioning between test platforms, Torontech’s application support team is available throughout the qualification process.
QC Consistent Instrument Performance & Traceability	Torontech instruments are manufactured under quality management systems aligned with ISO 9001 standards. Each unit ships with factory calibration documentation and traceable reference material certifications, ensuring your chemical analysis results are defensible in customer audits, third-party certifications, and regulatory inspections.
OEM Cost-Effective Performance Without Compromise	Torontech delivers laboratory-grade analytical performance at pricing that reflects a direct-to-customer supply model — providing world wide petroleum laboratories access to instruments that meet or exceed the specifications of major multinational brands, at a significantly lower total cost of ownership including instrument acquisition, consumables, and service.

Torontech’s Chemical Properties & Titration Analyzers are purpose-built for the petroleum and petrochemical industry. From crude oil intake to finished product release, our instruments deliver the analytical precision, method compliance, and operational reliability that modern petroleum laboratories demand — backed by North American application expertise and support.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

