



HOSPITAL BASSINET - INFANT BED

Mobile hospital bassinets and infant beds for newborn care, built to ISO 13485.

A hospital bassinet gives maternity wards, nurseries, and obstetrics units a safe, mobile bed for healthy newborns, keeping the infant in clear view and within reach of the care team. The ToronCare™ range spans four models, from a practical everyday unit to a premium adjustable infant bed, all built in line with ISO 13485

Each model pairs a clear tub with a wheeled, lockable frame, so staff can watch the baby, adjust position for feeding or reflux relief, and move the unit between rooms. The line ranges across frame materials, tilt systems, and height adjustment, which lets a facility match the bassinet to its ward, budget, and care routine.

Applications

Hospital bassinets support newborn care across the maternity and obstetrics pathway.

- Maternity and Postpartum Wards: keep the newborn at the mother's bedside for feeding, observation, and recovery.
- Newborn Nurseries: hold infants in clear-sided, mobile beds for grouped care and routine monitoring.
- Labor and Delivery Suites: receive the baby after birth and move with the mother to the recovery area.
- Obstetrics and Gynecology Units: position the infant for examination, diapering, and bedside care.
- Rooming-In and Mother-Baby Care: support family-centered care by keeping the infant beside the parent.
- Intra-Hospital Infant Transport: roll the newborn safely between wards, imaging, and treatment rooms on braked casters.

Hospital Bassinet Range

The ToronCare™ range covers four bassinets, each suited to a different ward, budget, and level of adjustability. Compare the models below.

Model	Frame & Basin	Tilt & Height	Best For
ToronCare™ 6200	Epoxy-coated steel, ABS tub	Gas-spring tilt 0-12°, fixed height	The everyday choice for routine newborn care and observation, with a wire storage basket and a compact footprint for tight wards.
ToronCare™ 6400	Stainless steel, ABS tub	Mechanical tilt 0-12°, fixed height	A stainless build for hygiene-focused, high-traffic wards; ships with a 4 cm leather-covered mattress and a stainless lower shelf.
ToronCare™ 6600	Aluminum column, food-grade PS tub	Gas-spring tilt 0-12°, foot-pedal height 910-1210 mm	A height-adjustable, column-mounted bassinet that raises and lowers hands-free for ergonomic access, with a push handle for transport.
ToronCare™ 6800	Welded steel (carbon + 304 stainless), PMMA acrylic basin	Two-way tilt $\pm 11^\circ$, hand-wrench height 840-970 mm	The premium model with two-way tilt and full positioning control, a scratch-resistant acrylic basin, and an under-frame protective grille.

How to Choose the Right Hospital Bassinet


- For routine newborn care on a budget: choose the ToronCare™ 6200, a practical epoxy-coated bassinet for everyday observation and transport.
- For hygiene-focused or high-traffic wards: choose the ToronCare™ 6400, a stainless steel model that stands up to constant cleaning and includes a mattress.
- For ergonomic, hands-free height adjustment: choose the ToronCare™ 6600, which raises and lowers by foot pedal to suit seated and standing caregivers.
- For full positioning control in obstetrics: choose the ToronCare™ 6800, with two-way tilt, adjustable height, and a premium acrylic basin.

Why Choose Torontech for Your Hospital Bassinet

Torontech North America supplies testing and measurement equipment to laboratories, quality teams, and manufacturers, backed by responsive technical support and calibration service. The hospital bassinet range is built for reliable, standards-driven results, so you can source confidently from a single partner. Our team helps you match the right configuration to your samples, applicable standards, and throughput, then supports you through installation and beyond.

Request a Quote

Tell us your ward type, the level of adjustment you need, and your volume, and we will recommend the right ToronCare™ hospital bassinet. [Contact the Torontech team](#) for pricing and lead times.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

