



MEDICAL CART - MEDICAL TROLLEY

Mobile medical carts for records, medication, emergency response, and clinical supply handling.

A medical trolley is a mobile cart that brings records, medications, supplies, or emergency equipment directly to the bedside. Hospitals, clinics, and long-term care facilities rely on these carts to keep clinical staff at the point of care rather than returning to a central room for every task.

This Torontech range spans medical record carts, medication carts, anesthesia and crash carts, utility and treatment trolleys, and soiled linen carts. Several models are built under ISO 13485 and ISO 9001 quality systems, so you can match the right cart to each ward while sourcing from one supplier.

Applications

Medical carts support nearly every clinical and support function in a hospital, from documentation and drug rounds to emergency response and waste removal. The right cart depends on the department and the task at hand.

- Inpatient Wards and Nursing Stations: hold patient charts, medications, and supplies at the bedside during rounds and shift handovers.
- Operating Rooms and Recovery: keep anesthesia drugs, surgical records, and resuscitation gear organized beside the surgical team.
- Emergency Departments and ICUs: stage crash carts and medication carts for fast access when a patient deteriorates.
- Outpatient Clinics and Day Surgery: move medications, files, and consumables between treatment bays for back-to-back caseloads.
- Treatment and Dressing Rooms: combine clean storage, a basin, and waste handling for wound care and minor procedures.
- Central Sterile and Supply Departments: transport instrument sets and consumables through cleaning, packing, and distribution.
- Housekeeping and Laundry: collect soiled linen and segregate waste streams during cleaning and changeover rounds.

Medical Cart Range

The range groups into seven cart families, each shaped around a clinical job. Use the groups below to find the cart type that fits your workflow, then compare the models within it.

Medical Record Carts

Mobile chart storage that holds patient files upright in divider bays and locks them between uses. These carts bring documentation to the bedside and support HIPAA-aligned handling at the unit level.

- ToronCAB™ 1100: single-bay patient chart cart with a lockable top drawer, a three-point cam-lock system, and four braking casters for ward rounds.
- ToronCAB™ 1150: double-bay version with two lockable columns and roll-down shutter doors for higher chart capacity or two staff working at once.

Medication Carts

Lockable-drawer carts for safe drug storage and bedside treatment, most with a single central lock and a sliding work surface. Capacity scales from a compact two-drawer unit to a five-drawer station, and each one keeps medications behind a lock between rounds.

- ToronCAB™ 1200: compact medication cart with two lockable drawers, an IV pole holder, a sliding table, and dual waste bins.
- ToronCAB™ 1230: standard cart adding a medicine cassette holder, a needle disposal container, and a disinfectant holder for safe sharps handling.
- ToronCAB™ 1280: locking medication cart with three drawers, an enclosed lower cabinet, and a central locking system for larger drug inventories.
- ToronCAB™ 1290: high-capacity cart with five drawers, an integrated IV pole, and a file cassette so documentation travels with the medications.

Anesthesia Cart

A perioperative cart that classifies drugs by color for fast, low-error selection during induction.

- ToronCAB™ 1300: five color-coded drawers split into sixteen adjustable compartments, with an upper medicine shelf, a side power strip, and 125 mm silent casters, built under ISO 13485.

Crash Carts

Emergency carts, also called crash carts or code carts, that put resuscitation drugs, a defibrillator shelf, oxygen, and a CPR board on one mobile platform. Breakaway seals show at a glance whether a hospital crash cart is sealed and fully stocked.

- ToronCAB™ 1420: compact crash cart with five color-coded drawers, a defibrillator shelf, a CPR board, and an oxygen cylinder holder.
- ToronCAB™ 1400: larger crash cart trolley that adds a file cassette for code records and more drawer and work space.
- ToronCAB™ 1460: premium HDPE crash cart with a lockable tabletop medicine box, an adjustable defibrillator tray, and a height-adjustable IV pole, built under ISO 13485.

Utility Carts

Open two-tier transport carts for instruments, supplies, and equipment, built in ABS or full stainless steel for frequently cleaned areas.

- ToronCAB™ 1500: two ABS decks on an epoxy-coated stainless steel frame, with guardrails and a 100 kg load capacity.
- ToronCAB™ 1530: full stainless steel utility cart with a side push handle and corner bumpers for sterile and high-traffic areas.
- ToronCAB™ 1550: stainless instrument cart that pairs two railed decks with an enclosed drawer for small tools and supplies.

Treatment and Dressing Carts

Self-contained carts for dressing changes and wound care, combining clean storage, a basin, and waste handling so a full procedure happens at one cart. Wipeable surfaces support infection-control routines between patients.

- ToronCAB™ 1600: ABS treatment trolley with a drawer, a stainless steel basin, and a waste bin under the decks.
- ToronCAB™ 1630: full stainless version of the same layout with a side push handle for sterile areas.
- ToronCAB™ 1650: wide three-tier stainless cart with a removable fabric waste bag for multi-patient rounds.

Soiled Linen Carts

Bagged carts that collect used linen and waste during bed changes and cleaning rounds, with removable bags for separate washing.

- ToronCAB™ 1700: single covered fabric bag in a stainless steel frame, with push handles at each end.
- ToronCAB™ 1730: dual-bag sorting cart with color-coded bags and soft-close lids for segregating linen and waste at the point of collection.

How to Choose the Right Medical Cart or Trolley

- For patient charts and records: choose a record cart, the single-bay ToronCAB™ 1100 or the double-bay 1300, with lockable drawers and divider bays.
- For medication rounds: choose a medication cart sized to your round, from the compact 2100 to the five-drawer 2700, each with central locking.
- For the operating room: choose the anesthesia cart (3500) with color-coded drawers and on-board power for monitors.
- For emergency and code response: choose a crash cart (4200 compact, 4400 larger, or 4600 premium) with a defibrillator shelf and breakaway seals.
- For general transport: choose a utility cart (5100 ABS or 5300 and 5500 stainless), with stainless steel preferred for sterile areas.


- For dressing and wound care: choose a treatment cart (6100 or 6300) with a basin and waste bin, or the wide 6500 for multi-patient rounds.
- For soiled linen and waste: choose a linen cart, the single-bag 7100 or the dual-bag sorting 7300.
- For sterile, frequently disinfected areas: choose full stainless steel construction over ABS wherever the cart offers it.

Why Choose Torontech for Your Medical Carts

Torontech supplies medical carts and clinical equipment to hospitals, clinics, and care facilities across North America and worldwide, backed by responsive technical support and service. The medical cart range is built for reliable, standards-driven daily use, so you can standardize on one supplier as your facility's needs grow. Our team helps you match the right configuration to your wards, workflows, and storage needs, then supports you through delivery and beyond.

Request a Quote

Tell us your department, the storage and locking each cart needs, and how many units you are outfitting, and we will recommend the right medical carts for your facility. [Contact the Torontech team](#) for pricing and lead times.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

