



MEDICAL CART (COMPACT)

Medical Cart (Compact) - ToronCAB™ 1200

ToronCAB™ 1200 compact medication cart organizes drugs, supplies, and infusion equipment for safe delivery at the bedside. Two lockable drawers, controlled by a single central lock, keep medications secured and separated from general supplies. A built-in IV pole holder supports infusion bags during treatment, while a sliding work surface gives staff room to prepare doses at the point of care. Four braking swivel casters and a push handle move the cart smoothly between rooms and hold it steady during use.

Nurses and clinical staff use the Torontech ToronCAB™ 1200 for medication rounds, treatment setups, and routine bedside care across hospital units. Two side-mounted waste bins let staff segregate sharps and general waste right at the cart, supporting safe disposal during procedures. The molded work surface and removable PVC top mat wipe down quickly between patients, which supports infection-control routines. Anti-bumper guards on all four corners protect walls, doorframes, and the cart through frequent moves.

For North American facilities, the central-locking drawers help keep medications secured between rounds, supporting controlled-access handling of drugs at the unit level. The compact L750 by W560 mm footprint moves through standard doorways and fits in tight treatment rooms and beside beds. Designed for daily clinical use, the ToronCAB™ 1200 brings medications, infusion support, a work surface, and waste disposal together on one mobile platform, so staff complete more of each task without leaving the patient.

APPLICATIONS

Medical Cart (Compact) - ToronCAB™ 1200 Applications

The ToronCAB™ 1200 supports medication delivery and bedside treatment across departments that administer drugs and fluids. Its drawers, IV support, and work surface fit a range of clinical settings.

- Hospital Inpatient Wards: carries medications and supplies for scheduled drug rounds, so nurses dispense at the bedside without repeated trips to the medication room.
- Treatment and Procedure Rooms: provides a mobile work surface and locked drug storage for dressing changes, injections, and minor procedures.
- Infusion and IV Therapy: holds infusion bags on the IV pole and keeps related supplies within reach during fluid and medication administration.
- Outpatient Clinics and Day Units: moves medications and consumables between treatment bays for back-to-back patient appointments.
- Intensive Care and High-Dependency Units: keeps frequently used drugs, fluids, and disposables organized and mobile in high-activity areas.
- Long-Term Care and Rehabilitation Facilities: supports medication rounds and routine treatments with secure storage in shared environments.
- Emergency and Recovery Areas: positions drugs, supplies, and a work surface near staff for fast medication preparation during active care.



FEATURES

Medical Cart (Compact) - ToronCAB™ 1200 Key Features

The ToronCAB™ 1200 combines secure medication storage, infusion support, and a mobile work surface in a compact medication cart built for daily clinical use.

- Two Lockable Drawers: drawers store drugs and supplies and stay organized at working height for quick access during rounds.
- Single Central Lock: one lock secures both drawers at once, so staff close down medications in a single action when leaving the cart.
- IV Pole Holder: an integrated pole supports infusion bags during treatment, keeping fluids ready and within reach.
- Sliding Work Surface: a pull-out table extends for dose preparation and retracts to keep the cart compact in transit.
- Dual Waste Bins: two side-mounted bins separate sharps and general waste at the point of use for safe disposal.
- Push Handle: a dedicated handle gives staff a firm grip and control when steering a loaded cart.
- Four Individual-Brake Casters: each swivel caster locks on its own for a stable base during use and free movement between rooms.
- Corner Anti-Bumper Guards: bumpers on all four corners cushion contact with walls and doorframes, protecting the facility and the cart.
- Removable PVC Top Mat: the protective mat lifts off for cleaning and shields the work surface from spills and scratches.

THEORY & METHOD

Theory and Method

The ToronCAB™ 1200 organizes medication delivery around secure, accessible drawer storage rather than open shelving. Two drawers hold drugs and supplies, and a single central lock secures both in one action, so staff close down the cart quickly when stepping away. Keeping medications behind a lock supports controlled access and reduces the chance of unattended drugs being disturbed. The drawers slide on the cart body at working height, so staff reach contents without bending or searching.

Treatment tasks need space and supporting equipment, which the cart provides through a sliding work surface, an IV pole holder, and waste bins. The sliding table extends to give staff a stable area to prepare doses, then retracts to keep the cart compact in transit. The IV pole holder secures infusion bags above the work area, so fluids hang ready during administration. Two waste bins let staff drop sharps and general waste into separate containers, which keeps disposal safe and immediate.

Mobility comes from four swivel casters, each with its own brake, guided by a dedicated push handle. The handle gives staff a firm grip and control when moving a loaded cart through corridors and doorways. Locking the four casters holds the cart steady while staff prepare or administer medications, and releasing them returns free movement in any direction. Corner anti-bumpers absorb knocks during transport, so the cart and the building avoid damage over years of daily rounds.

TECHNICAL SPECIFICATIONS


Medical Cart (Compact) - ToronCAB™ 1200 Technical Specifications

Parameter	Specification
External Dimensions	L750 × W560 × H1065 mm (±10 mm)
Medication Storage	Two lockable drawers with a single central lock
Work Surface	Molded top with removable PVC mat, plus a sliding side table
Infusion Support	Integrated IV pole holder

Parameter	Specification
Waste Management	Two side-mounted waste bins for sharps and general waste
Mobility	Four swivel casters with individual brakes

Standard Configuration

Item	Quantity	Description
Drawer	2 pcs	Lockable drawers for medications and supplies
Central Lock	1 pc	Single lock securing both drawers in one action
Waste Bin	2 pcs	Side-mounted bins for separating sharps and general waste
IV Pole Holder	1 pc	Holder supporting infusion bags during treatment
Sliding Working Table	1 pc	Pull-out work surface for dose preparation
Push Handle	1 pc	Handle for steering the loaded cart
Top Mat (PVC)	1 pc	Removable PVC mat protecting the top work surface
Anti-Bumper	4 pcs	Corner-mounted guards protecting walls, fixtures, and the cart
Individual Brake Caster	4 pcs	Swivel casters, each with a dedicated brake for stable positioning



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

