



TORONTECH

MEDICATION CART (SECURE LOCK MODEL)



Medication Cart (Secure Lock Model) - ToronCAB™ 1280

ToronCAB™ 1280 is a locking medication cart built for secure drug storage and bedside treatment on hospital wards. Three drawers and a lower cabinet store medications, fluids, and larger supplies, and a central locking system secures the drawers from a single point. A sliding work surface gives staff room to prepare doses at the point of care. Four braking swivel casters move the cart between rooms and hold it firmly in place during use.

Nurses and clinical staff use the Torontech ToronCAB™ 1280 for medication rounds and treatments where drug security and storage capacity both matter. The combination of three drawers and an enclosed cabinet holds a larger inventory than a basic trolley, so a single cart covers more of a round. A needle disposal container and two waste bins let staff segregate sharps and general waste at the bedside. Molded surfaces wipe down between patients, which supports infection-control routines on busy units.

For North American facilities, central locking supports controlled-access handling of medications, helping staff secure every drawer in one action before stepping away. The low 890 mm height and stable, wide base suit bedside work and rolling between rooms, while the cart still passes through standard doorways. Designed for daily clinical use, the ToronCAB™ 1280 brings secure multi-drawer storage, an enclosed cabinet, a work surface, and sharps disposal together on one platform, so staff manage medications safely throughout a shift.

APPLICATIONS

Medication Cart (Secure Lock Model) - ToronCAB™ 1280 Applications

The ToronCAB™ 1280 supports secure medication delivery and treatment across departments that handle a larger drug inventory or controlled medications. Its multi-drawer storage, locking cabinet, and sharps disposal fit a range of clinical settings.

- Hospital Inpatient Wards: carries a full round of medications in secured drawers, so nurses dispense at the bedside without repeated trips to the drug room.
- Medication and Drug-Storage Rooms: serves as a mobile, lockable extension of the pharmacy stock, keeping controlled and routine drugs organized and access-restricted.
- Treatment and Procedure Rooms: provides locked storage, a work surface, and sharps disposal for injections, dressing changes, and minor procedures.
- Intensive Care and High-Dependency Units: keeps a larger inventory of drugs and disposables organized and mobile in high-activity areas.
- Outpatient Clinics and Day Units: moves medicines and consumables between treatment bays with capacity for back-to-back appointments.
- Long-Term Care and Rehabilitation Facilities: supports medication rounds with secure, higher-capacity storage in shared environments.
- Recovery and Observation Areas: positions medicines, a work surface, and waste handling near staff for fast care after procedures.



FEATURES

Medication Cart (Secure Lock Model) - ToronCAB™ 1280 Key Features

The ToronCAB™ 1280 pairs a central locking system with multi-drawer and cabinet storage in a secure medication cart for daily clinical use.

- Central Locking System: one lock secures all three drawers at once, so staff close down the entire cart in a single action when leaving it.
- Three Storage Drawers: drawers separate medications by type or round and keep frequently used drugs at working height.
- Enclosed Lower Cabinet: a two-door cabinet adds covered storage for larger supplies behind a secured front.
- Sliding Work Surface: a pull-out table extends for dose preparation and retracts to keep the cart compact in transit.
- Needle Disposal Container: an on-board sharps container lets staff dispose of used needles immediately and safely at the bedside.
- Dual Waste Bins: two side-mounted bins separate general and clinical waste at the point of use.
- Four Individual-Brake Casters: each swivel caster locks on its own for a stable base during use and free movement between rooms.
- Wide, Low-Profile Build: the longer, lower body adds storage capacity and stability while still passing through standard doorways.

THEORY & METHOD

Theory and Method

The ToronCAB™ 1280 builds medication security around a central locking system rather than separate locks on each drawer. One lock controls all three drawers at once, so a single turn secures the whole cart when staff step away, and a single turn reopens it for the next patient. Closing everything in one action removes the risk of a drawer being left unlocked during a busy round. The enclosed lower cabinet adds covered storage for larger items behind the same secured front.

Storage on this model splits across three drawers and a sliding work surface, so staff separate medications by type or round and still have space to prepare doses. The three drawers stack at working height, keeping frequently used drugs reachable without bending or crouching. The sliding table extends to give a stable preparation area, then retracts to keep the cart compact in transit. Splitting medications across separate drawers also helps staff find the right item quickly during a fast-moving round.

Mobility comes from four swivel casters, each with its own brake. Locking all four casters holds the cart steady while staff work at the drawers or cabinet, and releasing them gives free movement in any direction through wards and doorways. A needle container and two waste bins mount on the cart, so staff drop sharps and general waste into the right container at the moment of use. Immediate, on-board disposal keeps used needles off work surfaces during procedures.

TECHNICAL SPECIFICATIONS


Medication Cart (Secure Lock Model) - ToronCAB™ 1280 Technical Specifications

Parameter	Specification
External Dimensions	L815 × W540 × H890 mm (±10 mm)
Medication Storage	Three drawers plus an enclosed lower cabinet
Locking	Central lock with central locking system for the drawers

Parameter	Specification
Work Surface	Sliding side table
Waste and Sharps	Two waste bins plus a needle disposal container
Mobility	Four swivel casters with individual brakes

Standard Configuration

Item	Quantity	Description
Drawer	3 pcs	Drawers for medications and supplies
Central Lock	1 pc	Lock mechanism for securing the drawers
Central Locking System	1 set	Secures all drawers from a single point
Sliding Working Table	1 pc	Pull-out work surface for medication preparation
Needle Disposal Container	1 pc	On-board sharps container for used needles
Waste Bin	2 pcs	Side-mounted bins for general and clinical waste
Individual Brake Caster	4 pcs	Swivel casters, each with a dedicated brake for stable positioning



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

