



MEDICATION CART (LARGE)

Medication Cart (Large) - ToronCAB™ 1290

ToronCAB™ 1290 hospital medication cart combines high-capacity drug storage, infusion support, and patient documentation in one mobile unit. Five drawers hold medications and supplies, secured by a central lock, while an integrated IV pole supports infusion bags during treatment. A file cassette keeps patient charts or paperwork with the cart, and a sliding work surface gives staff room to prepare doses. Four braking swivel casters move the cart between rooms and hold it steady at the bedside.

Nurses and clinical staff use the Torontech ToronCAB™ 1290 as a complete bedside station for medication rounds, infusion setups, and treatments. Five drawers give the cart more capacity than a basic two-drawer trolley, so staff carry a full round of drugs and consumables without restocking mid-round. A needle disposal container and two waste bins let staff segregate sharps and general waste at the point of use. Molded surfaces wipe down between patients, which supports infection-control routines on busy wards.

For North American facilities, the central lock helps keep medications secured between rounds, supporting controlled-access handling of drugs at the unit level. The L820 by W505 mm footprint moves through standard doorways and parks neatly beside beds and at nursing stations. Designed for daily clinical use, the ToronCAB™ 1290 brings drug storage, infusion support, patient paperwork, and sharps disposal together on one platform, so staff complete medication, treatment, and documentation tasks without leaving the patient.

APPLICATIONS

Medication Cart (Large) - ToronCAB™ 1290 Applications

The ToronCAB™ 1290 supports medication delivery, infusion, and bedside documentation across departments that need a higher-capacity, all-in-one cart. Its five drawers, IV pole, and file cassette fit a range of clinical settings.

- Hospital Inpatient Wards: carries a full round of medications across five drawers, so nurses dispense at the bedside without returning to the drug room to restock.
- Infusion and IV Therapy: holds infusion bags on the IV pole and keeps related supplies and paperwork together during fluid and medication administration.
- Treatment and Procedure Rooms: provides locked drug storage, a work surface, and sharps disposal for injections, dressing changes, and minor procedures.
- Intensive Care and High-Dependency Units: keeps a large inventory of drugs, fluids, and disposables organized and mobile in high-activity areas.
- Outpatient Clinics and Day Units: moves medicines, consumables, and patient files between treatment bays for back-to-back appointments.
- Long-Term Care and Rehabilitation Facilities: supports medication rounds and treatments with higher-capacity storage in shared environments.
- Recovery and Observation Areas: positions drugs, infusion support, a work surface, and waste handling near staff for fast care after procedures.



FEATURES

Medication Cart (Large) - ToronCAB™ 1290 Key Features

The ToronCAB™ 1290 pairs high-capacity, five-drawer storage with infusion support and patient documentation in a comprehensive hospital medication cart.

- Five Storage Drawers: more drawers separate medications by type, route, or round and let one cart carry a full medication round.
- Central Lock: a single lock secures all drawers, so staff close down the cart in one action when leaving it.
- Integrated IV Pole: an on-board pole supports infusion bags during treatment, keeping fluids ready and within reach.
- File Cassette: a holder keeps patient charts or treatment records with the cart, so documentation travels with the medications.
- Sliding Work Surface: a pull-out table extends for dose preparation and retracts to keep the cart compact in transit.
- Needle Disposal Container: an on-board sharps container lets staff dispose of used needles immediately and safely at the bedside.
- Dual Waste Bins: two side-mounted bins separate general and clinical waste at the point of use.
- Four Individual-Brake Casters: each swivel caster locks on its own for a stable base during use and free movement between rooms.

THEORY & METHOD

Theory and Method

The ToronCAB™ 1290 organizes a large medication inventory across five separate drawers rather than a few deep compartments. Splitting storage this way lets staff group drugs by type, route, or round, so the right item is easy to find during a fast-moving shift. A single central lock secures all the drawers, so staff close down the cart in one action when stepping away. The drawers stack at working height, keeping frequently used medications reachable without bending.

Beyond storage, the cart supports the wider tasks of a medication round through an IV pole, a file cassette, and a sliding work surface. The IV pole holds infusion bags above the work area, so fluids hang ready while staff prepare related medications. The file cassette keeps patient charts or treatment records with the cart, so documentation travels alongside the drugs it relates to. The sliding table extends for dose preparation and retracts to keep the cart compact in transit.

Mobility comes from four swivel casters, each with its own brake. Locking all four casters holds the loaded cart steady while staff work at the drawers, and releasing them gives free movement in any direction through wards and doorways. A needle disposal container and two waste bins mount on the cart, so staff drop sharps and general waste into the right container at the moment of use. Keeping disposal on board removes the need to carry used needles across the room.

TECHNICAL SPECIFICATIONS

Medication Cart (Large) - ToronCAB™ 1290 Technical Specifications

Parameter	Specification
External Dimensions	L820 × W505 × H1000 mm (±10 mm)
Medication Storage	Five drawers secured by a central lock
Documentation	File cassette for patient charts or records
Infusion Support	Integrated IV pole
Work Surface	Sliding side table

Parameter	Specification
Waste and Sharps	Two waste bins plus a needle disposal container
Mobility	Four swivel casters with individual brakes

Standard Configuration

Item	Quantity	Description
Drawer	5 pcs	Drawers for medications and supplies
Central Lock	1 pc	Lock securing all drawers from one point
Sliding Working Table	1 pc	Pull-out work surface for medication preparation
I.V. Pole	1 pc	Pole supporting infusion bags during treatment
File Cassette	1 pc	Holder for patient charts or treatment records
Needle Disposal Container	1 pc	On-board sharps container for used needles
Waste Bin	2 pcs	Side-mounted bins for general and clinical waste
Individual Brake Caster	4 pcs	Swivel casters, each with a dedicated brake for stable positioning



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

