



TORONTECH

UTILITY CART (ABS FEATURE)



Utility Cart (ABS Feature) - ToronCAB™ 1500

ToronCAB™ 1500 is a two-tier medical utility cart for moving instruments, supplies, and equipment around wards, clinics, and labs. Two ABS (Anti-lock Braking System) plastic decks sit on an epoxy-coated stainless steel frame, giving staff organized storage on a surface that wipes clean between uses. Guardrails around each deck keep items from sliding off in transit. Four 100 mm casters, two with brakes, roll the cart smoothly and hold it steady when parked, while the frame carries loads up to 100 kg.

Hospitals and clinics choose the Torontech ToronCAB™ 1500 as a simple, durable workhorse for everyday transport tasks. The stainless steel frame resists corrosion in damp, frequently cleaned environments, and the molded ABS decks hold up to repeated wiping and disinfection. The cart ships in knock-down form, which lowers freight cost and lets staff assemble it on site. With open shelving and no fixed function, the cart adapts to whatever a department needs to move.

The compact 720 by 460 mm footprint moves easily through corridors and doorways and tucks into tight treatment or storage areas. Built for daily use, the ToronCAB™ 1500 gives any department a low-cost, easy-to-clean cart for organized, mobile storage.

APPLICATIONS

Utility Cart (ABS Feature) - ToronCAB™ 1500 Applications

The ToronCAB™ 1500 suits any department that needs to move instruments, supplies, or equipment between locations. Its open two-tier design adapts to a wide range of clinical and support settings.

- Treatment and Examination Rooms: carries instruments, dressings, and consumables to the bedside and between exam rooms.
- Wards and Nursing Stations: moves supplies, linens, or equipment during rounds and restocking.

- Operating Room Support: transports sterile packs and instruments between storage, prep, and theater areas.
- Laboratories: carries samples, glassware, and equipment between benches and storage.
- Central Sterile and Supply Departments: moves instrument sets and consumables through cleaning, packing, and distribution.
- Dental and Outpatient Clinics: keeps tools and materials organized and mobile during procedures.
- Pharmacies and Storerooms: relocates stock and supplies during picking and restocking.



FEATURES

Utility Cart (ABS Feature) - ToronCAB™ 1500 Key Features

The ToronCAB™ 1500 combines a durable stainless steel frame, wipeable ABS decks, and easy mobility in a simple two-tier utility cart.

- Two ABS Plastic Decks: two molded shelves give organized, wipeable storage on a compact footprint.
- Epoxy-Coated Stainless Steel Frame: a corrosion-resistant frame stands up to damp, frequently cleaned hospital environments.
- 100 kg Load Capacity: the frame carries up to 100 kg across the two decks for heavier transport tasks.
- Guardrails on Each Deck: rails around both shelves keep items from sliding off in transit.
- Four 100 mm Casters with Brakes: four casters roll smoothly over hospital floors, and two with brakes hold the cart steady when parked.
- Knock-Down Construction: the cart ships flat to lower freight cost and assembles on site without special tools.

THEORY & METHOD

Theory and Method


The ToronCAB™ 1500 keeps transport simple by separating loads across two open decks rather than enclosing them in drawers. Two tiers double the usable surface within a small footprint, so staff group items by deck and load or unload from any side. A guardrail runs around each deck to stop items sliding off as the cart rolls over thresholds and around corners. Open shelving also lets staff see and reach everything at a glance, without opening or searching compartments.

Durability comes from the materials and frame. An epoxy-coated stainless steel frame resists rust and supports loads up to 100 kg, while molded ABS decks give smooth, wipeable surfaces that stand up to cleaning. Four 100 mm casters carry the cart over hospital floors, and brakes on two of them lock it in place during loading or use. The frame assembles from a knock-down kit, which ships flat to cut freight cost and goes together on site without special tools.

TECHNICAL SPECIFICATIONS

Utility Cart (ABS Feature) - ToronCAB™ 1500 Technical Specifications

Parameter	Specification
External Dimensions	720 (L) × 460 (W) × 970 (H) mm
Frame	Epoxy-coated stainless steel
Decks	Two ABS engineering plastic decks
Guardrails	Stainless steel rail around each deck
Casters	Four × 100 mm; two with brakes
Load Capacity	100 kg
Construction	Knock-down (flat-pack)



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

