



TORONTECH



UTILITY CART (PREMIUM)

Utility Cart (Premium) - ToronCAB™ 1550

ToronCAB™ 1550 is a stainless steel instrument cart that pairs two open shelves with an enclosed drawer for organized instrument transport. Both decks carry guardrails to keep trays and instruments from sliding off, while a drawer under the lower deck stores smaller tools and supplies out of sight. Full stainless steel construction gives a hygienic surface that withstands repeated cleaning and disinfection. Four 100 mm casters, two with brakes, roll the cart smoothly and hold it steady when parked.

Hospitals and clinics choose the Torontech ToronCAB™ 1550 when a cart needs both open trays and a place to keep small items secure. Staff lay instruments or dressing packs on the decks and store gauze, blades, or fittings in the drawer, so everything for a procedure travels together. Stainless steel resists corrosion and tolerates strong disinfectants, which makes the cart a good fit for treatment rooms, clinics, and sterile areas. The compact 665 by 475 mm footprint moves easily through corridors and tight spaces.

Built for daily use, the ToronCAB™ 1550 gives treatment and procedure areas a durable, easy-to-clean instrument cart that keeps both open and enclosed storage on one mobile platform.

APPLICATIONS

Utility Cart (Premium) - ToronCAB™ 1550 Applications

The ToronCAB™ 1550 suits procedures and tasks that need both open shelves and a closed drawer for instruments and small supplies. Its stainless build fits clinical and sterile settings.

- Treatment and Procedure Rooms: holds instruments and dressing packs on the decks and keeps small tools in the drawer for minor procedures.

- Operating Room Support: moves sterile trays and instruments between storage, prep, and theater, with a drawer for smaller items.
- Dressing and Wound Care: keeps dressings and instruments on the decks and consumables in the drawer at the bedside.
- Dental and Outpatient Clinics: organizes tools on the shelves and small materials in the drawer during procedures.
- Examination Rooms: carries instruments and supplies to the bedside and between rooms.
- Laboratories: moves samples and equipment on the decks and keeps small items secured in the drawer.
- Central Sterile and Supply Departments: transports instrument sets through cleaning, packing, and distribution.



FEATURES

Utility Cart (Premium) - ToronCAB™ 1550 Key Features

The ToronCAB™ 1550 combines open shelving, an enclosed drawer, and durable stainless construction in a compact instrument cart.

- Two Decks with Guardrails: two open stainless shelves with rails on each keep trays and instruments organized and stop them sliding off in transit.
- Drawer Under the Lower Deck: an enclosed drawer stores smaller tools, sharps, or supplies that staff want contained rather than loose on a shelf.
- Full Stainless Steel Construction: frame, decks, and drawer are stainless steel, giving a corrosion-resistant surface that tolerates strong disinfection.
- Four 100 mm Casters with Brakes: four casters roll smoothly over hospital floors, and two with brakes hold the cart steady when parked.
- Compact Footprint: the 665 by 475 mm base keeps the cart maneuverable in tight treatment rooms and corridors.

THEORY & METHOD

Theory and Method

The ToronCAB™ 1550 combines two kinds of storage on one frame so staff keep open and enclosed items together. Two open decks give quick, visible access to trays and instruments, with a guardrail around each to stop items sliding off during transport. A drawer under the lower deck adds enclosed space for smaller tools, sharps, or supplies that staff want contained rather than loose on a shelf. Splitting storage this way lets one cart carry everything a procedure needs in a logical layout.

Material and mobility complete the design. Full stainless steel construction gives a single cleanable surface with no plastic to crack or stain, so the cart holds up to strong disinfection in sterile and frequently cleaned areas. Four 100 mm casters carry the cart over hospital floors, and brakes on two of them lock it in place while staff work at the decks or drawer. The compact frame keeps the cart maneuverable in tight treatment rooms and crowded corridors.

TECHNICAL SPECIFICATIONS

Utility Cart (Premium) - ToronCAB™ 1550 Technical Specifications

Parameter	Specification
External Dimensions	665 (L) × 475 (W) × 990 (H) mm (±10 mm)
Construction	Full stainless steel (frame, decks, and drawer)
Decks	Two stainless steel decks with guardrails
Drawer	One drawer under the lower deck
Casters	Four × 100 mm; two with brakes



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

