

## LINEN CART (DOUBLE)



## **Linen Cart (Double) - ToronCAB™ 1730**

ToronCAB™ 1730 is a soiled linen sorting cart that holds two color-coded bags, so staff separate linen and waste streams at the point of collection. Each bag sits in a stainless steel frame and takes a disposable plastic liner for clean, sanitary disposal. A hard engineering-plastic lid covers each bag, and a buffer device lets the lid open and close quietly without slamming. Four 100 mm casters, two with brakes, and a push handle move the cart between rooms and hold it steady.

Hospitals and clinics use the Torontech ToronCAB™ 1730 to segregate soiled items as they collect them, which supports infection-control and laundry-sorting routines. Color-coded bags let staff route clean linen, soiled linen, or clinical waste into the right bag at the bedside, so sorting happens once instead of later. Disposable plastic liners keep each bag sanitary and make emptying quick, while the covered lids contain odor and shield the contents during transport. The stainless steel frame resists corrosion and wipes down for regular cleaning.

The 905 by 455 mm footprint keeps the dual-bag cart compact enough for corridors and ward bays, and the engineering-plastic lids hold up to daily use. Built for steady collection rounds, the ToronCAB™ 1730 brings sorting, hygiene, and mobility together in one stainless soiled linen cart.

### **APPLICATIONS**

---

## **Linen Cart (Double) - ToronCAB™ 1730 Applications**

The ToronCAB™ 1730 suits settings that separate linen and waste at the point of collection rather than sorting it afterward. Its two color-coded bags fit a range of hospital and facility routines.

- Hospital Wards: separates clean and soiled linen, or linen and waste, during bed changes and rounds.
- Infection Control and Waste Segregation: routes contaminated and general waste into separate bags at the bedside.
- Laundry Sorting: divides linen by type or soil level so the laundry receives pre-sorted loads.

- Operating Rooms and Recovery: separates used drapes, gowns, and waste during turnover between cases.
- Housekeeping and Cleaning Rounds: keeps recyclable, general, and soiled items apart while staff clean.
- Care Homes and Rehabilitation Facilities: supports routine linen changes with sorted collection.



## FEATURES

---

### Linen Cart (Double) - ToronCAB™ 1730 Key Features

The ToronCAB™ 1730 combines dual-bag sorting, soft-close lids, and a durable stainless frame in a single soiled linen cart.

- Two Color-Coded Sorting Bags: two fabric bags let staff separate linen and waste streams at the point of collection.
- Disposable Plastic Liners: each bag takes a plastic liner for sanitary disposal and quick, no-contact emptying.
- Engineering-Plastic Lids with Soft Close: a hard lid covers each bag and closes on a buffer device, shutting quietly to contain odor and contents.
- Stainless Steel Frame: a corrosion-resistant frame holds both loaded bags and wipes clean after each use.
- Push Handle: a handle gives staff a firm grip and control when steering the loaded cart.
- Four 100 mm Casters with Brakes: four casters roll smoothly over hospital floors, and two with brakes hold the cart steady when parked.

## THEORY & METHOD

---

### Theory and Method

The ToronCAB™ 1730 turns collection and sorting into a single step by giving staff two bags instead of one. Each bag carries a different stream, so a worker drops clean linen, soiled linen, or waste into the matching bag, and nothing has to be re-sorted later. Color coding makes the choice obvious at a glance, which cuts mistakes during a busy round. A disposable plastic liner lines each bag, so emptying means lifting out a sealed liner rather than handling soiled contents directly.

Construction adds hygiene and durability to the sorting design. A hard engineering-plastic lid covers each bag and closes on a buffer device, so it shuts quietly and stays down to contain odor and contents. The stainless steel frame holds both loaded

bags, resists corrosion, and wipes clean after each use. Four 100 mm casters carry the cart over hospital floors with brakes on two for stability, and a push handle gives staff a firm grip when steering the loaded cart.

TECHNICAL SPECIFICATIONS

Linen Cart (Double) - ToronCAB™ 1730 Technical Specifications

| Parameter           | Specification                                                     |
|---------------------|-------------------------------------------------------------------|
| External Dimensions | 905 (L) × 455 (W) × 905 (H) mm (±10 mm)                           |
| Frame               | Stainless steel                                                   |
| Bags                | Two color-coded fabric bags (with disposable plastic liners)      |
| Lids                | Engineering-plastic top cover on each bag, with soft-close buffer |
| Handle              | Push handle                                                       |
| Casters             | Four × 100 mm; two with brakes                                    |
| Carton Size         | 940 × 630 × 920 mm (1 pc/box)                                     |
| Gross Weight        | 19 kg                                                             |
| Shipping Volume     | 0.54 m³                                                           |



**Torontech Inc.**

251 Consumers Road, Suite 1200  
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652  
Email: [info@torontech.com](mailto:info@torontech.com)

**Torontech, LLC**

601 Brickell Key Drive, Suite 700  
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919  
Email: [sales@torontech.com](mailto:sales@torontech.com)

**Torontech ME FZE**

Meydan Grandstand, 6th Floor,  
Meydan Road, Nad Al Sheba,  
Dubai-261440, UAE 

Phone: +971 488 19252  
Email: [me@torontech.com](mailto:me@torontech.com)

